

Datasheet: STAR128

BATCH NUMBER 163465

Description:	GOAT F(ab') ₂ ANTI HUMAN KAPPA LIGHT CHAIN
Specificity:	KAPPA LIGHT CHAIN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Human

Product Form F(ab)₂ fragments of purified IgG - liquid

Antiserum Preparation Antisera to human kappa chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on pooled human Igs with kappa light chains covalently linked to agarose. F(ab)₂ fragments were prepared by pepsin digestion.

Buffer Solution Borate buffered saline

Preservative Stabilisers None present

Approx. Protein F(ab')₂ concentration 0.5 mg/ml

Concentrations

Immunogen Human kappa light chains.

External Database Links

UniProt:

[P01834](#) [Related reagents](#)

Entrez Gene:

[3514](#) IGKC [Related reagents](#)

RRID AB_1102712

Specificity **Goat F(ab')₂ anti Human Kappa Light Chain antibody** recognizes human kappa light chains. This antibody has been cross absorbed with pooled human myeloma proteins with lambda light chains.

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10123 available at:
<https://www.bio-rad-antibodies.com/SDS/STAR128>
10123

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M401861:220718'

Printed on 13 Aug 2023
