

Datasheet: STAR128PE

BATCH NUMBER 166655

| | |
|----------------------|---|
| Description: | GOAT F(ab') ₂ ANTI HUMAN KAPPA LIGHT CHAIN:RPE |
| Specificity: | KAPPA LIGHT CHAIN |
| Format: | RPE |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 0.5 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry | ■ | | | 1/20 - 1/30 |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

| | | | |
|----------------|--|---------------------|-------------------|
| Target Species | Human | | |
| Product Form | F(ab)2 fragments of purified IgG conjugated to R. Phycoerythrin - liquid | | |
| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) |
| | RPE 488nm laser | 496 | 578 |

Antiserum Preparation Antisera to human kappa chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography. F(ab)₂ fragments were prepared by pepsin digestion.

Buffer Solution Phosphate buffered saline

Preservative <0.1% Sodium Azide (NaN₃)
Stabilisers Stabilizing agent (sucrose)

Approx. Protein Concentrations F(ab')₂ concentration 0.25 mg/ml

Immunogen Human kappa light chains.

External Database

Links

UniProt:

[P01834](#)

[Related reagents](#)

Entrez Gene:

[3514](#)

IGKC

[Related reagents](#)

RRID

AB_1102709

Specificity

Goat F(ab')₂ anti Human Kappa Light Chain antibody recognizes human kappa light chains. This antibody has been cross absorbed with pooled human myeloma proteins with lambda light chains.

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul

Storage

Store at +4°C.

DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

6 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10045 available at:
<https://www.bio-rad-antibodies.com/SDS/STAR128PE>
10045

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M421242:230706'

Printed on 01 Mar 2024