

## Datasheet: STAR129

**BATCH NUMBER 167432**

|                      |                                    |
|----------------------|------------------------------------|
| <b>Description:</b>  | GOAT ANTI HUMAN LAMBDA LIGHT CHAIN |
| <b>Specificity:</b>  | LAMBDA LIGHT CHAIN                 |
| <b>Format:</b>       | Purified                           |
| <b>Product Type:</b> | Polyclonal Antibody                |
| <b>Isotype:</b>      | Polyclonal IgG                     |
| <b>Quantity:</b>     | 1 mg                               |

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                            | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry             |     |    | ▪              |                    |
| Immunohistology - Frozen   |     |    | ▪              |                    |
| Immunohistology - Paraffin |     |    | ▪              |                    |
| ELISA                      | ▪   |    |                |                    |
| Immunoprecipitation        |     |    | ▪              |                    |
| Western Blotting           |     |    | ▪              |                    |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**Target Species** Human

**Product Form** Purified IgG - liquid

**Antiserum Preparation** Antisera to human lambda chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on pooled human igs with lambda light chains covalently linked to agarose.

**Buffer Solution** Borate buffered saline

**Preservative Stabilisers** None present

**Approx. Protein** IgG concentration 1.0 mg/ml

## Concentrations

---

**Immunogen** Human lambda light chains.

---

## External Database Links

### UniProt:

[P0CG04](#) [Related reagents](#)

### Entrez Gene:

[28815](#) IGLV2-14 [Related reagents](#)

---

**RRID** AB\_1102720

---

**Specificity** **Goat anti Human Lambda Light Chain antibody** recognizes human lambda light chains. This antibody has been cross absorbed with pooled human myeloma proteins bearing kappa light chains.

---

**References**

1. Ben-horin, S. *et al.* (2016) Cross-immunogenicity: antibodies to infliximab in Remicade-treated patients with IBD similarly recognise the biosimilar Remsima. [Gut. 65 \(7\): 1132-8.](#)
2. Magro, F. *et al.* (2018) The performance of Remicade<sup>®</sup>- optimized quantification assays in the assessment of Flixabi<sup>®</sup> levels. [Therap Adv Gastroenterol. 11: 1756284818796956.](#)

---

**Storage** Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

---

**Health And Safety Information** Material Safety Datasheet documentation #10123 available at: <https://www.bio-rad-antibodies.com/SDS/STAR129>  
10123

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M421243:230706'

Printed on 01 Mar 2024

