

Datasheet: STAR130

BATCH NUMBER 171480

Description:	GOAT F(ab') ₂ ANTI HUMAN LAMBDA LIGHT CHAIN
Specificity:	LAMBDA LIGHT CHAIN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Human

Product Form F(ab')₂ fragments of purified IgG - liquid

Antiserum Preparation Antisera to human lambda chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography. F(ab')₂ fragments were prepared by pepsin digestion.

Buffer Solution Borate buffered saline

Preservative Stabilisers None present

Approx. Protein Concentrations F(ab')₂ concentration 0.5mg/ml

Immunogen Human lambda light chains.

External Database

Links

UniProt:

[P0CG04](#) [Related reagents](#)

Entrez Gene:

[28815](#) IGLV2-14 [Related reagents](#)

RRID

AB_1102722

Specificity

Goat F(ab')₂ anti Human Lambda Light Chain antibody recognizes human lambda light chains. This antibody has been cross absorbed against pooled human myeloma proteins with kappa light chains.

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information

Material Safety Datasheet documentation #10123 available at:
<https://www.bio-rad-antibodies.com/SDS/STAR130>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M428674:240301'

Printed on 29 Jan 2026
