

Datasheet: STAR132

Description:	GOAT ANTI MOUSE IgG1
Specificity:	IgG1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

RRID AB_1102653

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Mouse

Species Cross Reactivity Does not react with:Human

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to mouse IgG1 were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution Borate buffered saline

Preservative Stabilisers 0.09% Sodium Azide

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen Mouse IgG1 paraproteins.

External Database Links	<p>UniProt:</p> <p>P01869 Related reagents</p> <p>P01868 Related reagents</p> <p>Entrez Gene:</p> <p>16017 Ighg1 Related reagents</p> <p>16017 Ighg1 Related reagents</p>
Synonyms	Igh-4
Specificity	<p>Goat anti Mouse IgG1 antibody recognizes Mouse IgG1. This antibody has been cross absorbed against mouse IgM, IgG2a, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG1 antibody shows minimal cross-reactivity with human immunoglobulins.</p>
References	<ol style="list-style-type: none"> 1. Croft, N.P. <i>et al.</i> (2009) Stage-specific inhibition of MHC class I presentation by the Epstein-Barr virus BNLF2a protein during virus lytic cycle. PLoS Pathog. 5(6): e1000490. 2. Zuo, J. <i>et al.</i> (2011) The Epstein-Barr virus-encoded BILF1 protein modulates immune recognition of endogenously processed antigen by targeting MHC class I molecules trafficking on both the exocytic and endocytic pathways. J Virol. 85: 1604-14. 3. Knipping, K. <i>et al.</i> (2011) A gastrointestinal rotavirus infection mouse model for immune modulation studies. Virology 453: 109. 4. Young, D. <i>et al.</i> (2012) Soy-derived di- and tripeptides alleviate colon and ileum inflammation in pigs with dextran sodium sulfate-induced colitis. J Nutr. 142 (2): 363-8. 5. Bagai, U. and Pawar, A. (2013) A blood stage fraction of <i>Plasmodium berghei</i> induces protective and long lasting immune response in BALB/c mice. Parasitol Int. 62: 329-36. 6. Anda, S. <i>et al.</i> (2014) Cell-cycle analyses using thymidine analogues in fission yeast. PLoS One. 9 (2): e88629. 7. Kamat, M.M. <i>et al.</i> (2016) Changes in Myeloid Lineage Cells in the Uterus and Peripheral Blood of Dairy Heifers During Early Pregnancy. Biol Reprod. Aug 10. pii: biolreprod.116.141069. [Epub ahead of print] 8. Ramanathan, R. <i>et al.</i> (2015) Transplantation of human stem cell-derived hepatocytes in an animal model of acute liver failure. Surgery. 158 (2): 349-59. 9. Hwang, S.R. <i>et al.</i> (2015) Altered expression levels of neurodevelopmental proteins in fetal brains of BTBR T+tf/J mice with autism-like behavioral characteristics. J Toxicol Environ Health A. 78 (8): 516-23. 10. Zhao, Z. <i>et al.</i> (2015) Multiple B-cell epitope vaccine induces a Staphylococcus enterotoxin B-specific IgG1 protective response against MRSA infection. Sci Rep. 5: 12371.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10077 available at: 10077: https://www.bio-rad-antibodies.com/uploads/MSDS/10077.pdf

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M339659:181219'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)