

## Datasheet: STAR132P

**BATCH NUMBER 167817**

<b>Description:</b>	GOAT ANTI MOUSE IgG1:HRP
<b>Specificity:</b>	IgG1
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA	▪			1/4000 - 1/8000
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

#### Target Species

Mouse

#### Species Cross Reactivity

Does not react with:Human

#### Product Form

Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

#### Antiserum Preparation

Antisera to murine IgG1 were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared by affinity chromatography.

#### Buffer Solution

Phosphate buffered saline

#### Preservative Stabilisers

0.01% Thiomersal  
50% Glycerol

#### Immunogen

Mouse IgG1 paraproteins.

---

**External Database****Links****UniProt:**[P01869](#)    [Related reagents](#)[P01868](#)    [Related reagents](#)**Entrez Gene:**[16017](#) Ighg1    [Related reagents](#)[16017](#) Ighg1    [Related reagents](#)

---

**Synonyms**

Igh-4

---

**RRID**

AB\_2124272

---

**Specificity**

**Goat anti Mouse IgG1 antibody** recognizes Mouse IgG1. This antibody has been cross absorbed against mouse IgM, IgG2a, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG1 antibody shows minimal cross-reactivity with human immunoglobulins.

Goat anti Mouse IgG1 antibody has been validated for use on the [Genesis Cell Isolation System with the CelSelect Slide™ technology](#).

---

**References**

1. Croft, N.P. *et al.* (2009) Stage-specific inhibition of MHC class I presentation by the Epstein-Barr virus BNLF2a protein during virus lytic cycle. [PLoS Pathog. 5\(6\): e1000490.](#)
2. Zuo, J. *et al.* (2011) The Epstein-Barr virus-encoded BILF1 protein modulates immune recognition of endogenously processed antigen by targeting MHC class I molecules trafficking on both the exocytic and endocytic pathways. [J Virol. 85: 1604-14.](#)
3. Knipping, K. *et al.* (2011) A gastrointestinal rotavirus infection mouse model for immune modulation studies. [Virology J. 8: 109.](#)
4. Young, D. *et al.* (2012) Soy-derived di- and tripeptides alleviate colon and ileum inflammation in pigs with dextran sodium sulfate-induced colitis. [J Nutr. 142 \(2\): 363-8.](#)
5. Bagai, U. and Pawar, A. (2013) A blood stage fraction of *Plasmodium berghei* induces protective and long lasting immune response in BALB/c mice. [Parasitol Int. 62: 329-36.](#)
6. Anda, S. *et al.* (2014) Cell-cycle analyses using thymidine analogues in fission yeast. [PLoS One. 9 \(2\): e88629.](#)
7. Kamat, M.M. *et al.* (2016) Changes in Myeloid Lineage Cells in the Uterus and Peripheral Blood of Dairy Heifers During Early Pregnancy. [Biol Reprod. 95 \(3\): 68.](#)
8. Ramanathan, R. *et al.* (2015) Transplantation of human stem cell-derived hepatocytes in an animal model of acute liver failure. [Surgery. 158 \(2\): 349-59.](#)
9. Hwang, S.R. *et al.* (2015) Altered expression levels of neurodevelopmental proteins in fetal brains of BTBR T+tf/J mice with autism-like behavioral characteristics. [J Toxicol Environ Health A. 78 \(8\): 516-23.](#)
10. Zhao, Z. *et al.* (2015) Multiple B-cell epitope vaccine induces a Staphylococcus enterotoxin B-specific IgG1 protective response against MRSA infection. [Sci Rep. 5: 12371.](#)
11. Gutiérrez-Miranda, B. *et al.* (2020) Oleacein Attenuates the Pathogenesis of Experimental Autoimmune Encephalomyelitis through Both Antioxidant and Anti-Inflammatory Effects. [Antioxidants \(Basel\). 9 \(11\): 1161.](#)

12. Gutierrez, B. *et al.* (2020) Oleanolic acid ameliorates intestinal alterations associated with EAE. [J Neuroinflammation. 17 \(1\): 363.](#)
13. Apóstolo, N. *et al.* (2020) Synapse type-specific proteomic dissection identifies IgSF8 as a hippocampal CA3 microcircuit organizer. [Nat Commun. 11 \(1\): 5171.](#)
14. Zhuang, X. *et al.* (2020) CAR T cells targeting tumor endothelial marker CLEC14A inhibit tumor growth. [JCI Insight. 5 \(19\): e138808.](#)
15. Sparks, A.M. *et al.* (2018) Natural Selection on Antihelminth Antibodies in a Wild Mammal Population. [Am Nat. 192 \(6\): 745-760.](#)
16. Kushwaha, V. *et al.* (2019) Troponin 1 of human filarial parasite *Brugia malayi*: cDNA cloning, expression, purification, and its immunoprophylactic potential. [Parasitol Res. 118 \(6\): 1849-63.](#)
17. Kim, Y.H. *et al.* (2019) Effects of ruminal pH on gene expression in the rumen epithelium, peripheral blood mononuclear cell subpopulations, and blood metabolites from Holstein calves during weaning transition. [J Vet Med Sci. 81 \(6\): 808-816.](#)
18. Jimbo, S. *et al.* (2019) Natural and inducible regulatory B cells are widely distributed in ovine lymphoid tissues. [Vet Immunol Immunopathol. 211: 44-48.](#)
19. Sparks, A.M. *et al.* (2019) The genetic architecture of helminth-specific immune responses in a wild population of Soay sheep (*Ovis aries*). [PLoS Genet. 15 \(11\): e1008461.](#)
20. Nedumpun, T. *et al.* (2019) Negative Immunomodulatory Effects of Type 2 Porcine Reproductive and Respiratory Syndrome Virus-Induced Interleukin-1 Receptor Antagonist on Porcine Innate and Adaptive Immune Functions. [Front Immunol. 10: 579.](#)
21. Verma, R. *et al.* (2018) *Leishmania donovani*. molecules recognized by sera of filaria infected host facilitate filarial infection. [Parasitol Res. 117 \(9\): 2901-12.](#)
22. Gatkowska, J. *et al.* (2019) The Impact of the Antigenic Composition of Chimeric Proteins on Their Immunoprotective Activity against Chronic Toxoplasmosis in Mice. [Vaccines \(Basel\). 7\(4\):154.](#)
23. Han, H.J. *et al.* (2022) Metal arsenic mediated enhancement of type-2 immunity in brains with altered locomotive activities in mice with autism-like behavioral characteristics. [Toxicol Res. 38 \(1\): 27-33.](#)
24. Anda, S. *et al.* (2021) Cosegregation of asymmetric features during cell division. [Open Biol. 11 \(8\): 210116.](#)
25. Minaei, S. *et al.* (2018) Propranolol efficacy as a novel adjuvant for immunization against *Toxoplasma gondii*. tachyzoites. [Exp Parasitol. 194: 60-66.](#)
26. Ngoepe, T.A. *et al.* (2023) Th1 and Th2 epitopes of Cowdria polymorphic gene 1 of *Ehrlichia ruminantium*. [Onderstepoort J Vet Res. 90 \(1\): e1-e15.](#)
27. Jo, J. *et al.* (2023) Immunodysregulatory potentials of polyethylene or polytetrafluorethylene microplastics to mice subacutely exposed via intragastric intubation. [Toxicol Res. 39 \(3\): 419-27.](#)
28. Sessevmez, M. *et al.* (2023) Induction of humoral and cell-mediated immunity in mice by chitosan-curdlan composite nanoparticles administered intranasally and subcutaneously [J Drug Deliv Sci Technol. 86: 104704.](#)
29. Veisi, R. *et al.* (2024) Expression and purification of TolC as a recombinant protein vaccine against *Shigella flexneri* and evaluation of immunogenic response in mice. [Microb Pathog. : 106539 \[Epub ahead of print\]](#)
30. Faber, E. *et al.* (2024) Identification of T cell and linear B cell epitopes on African horse sickness virus serotype 4 proteins VP1-1, VP2, VP4, VP7 and NS3. [Vaccine. 42 \(2\):](#)

[136-145.](#)

31. Pretorius, A. *et al.* (2024) Screening for immune biomarkers associated with infection or protection against *Ehrlichia ruminantium* by RNA-sequencing analysis. [Microb Pathog. : 106588.](#)

---

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10097 available at: <https://www.bio-rad-antibodies.com/SDS/STAR132P>  
10097

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M421255:230706'

Printed on 09 Mar 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)