

Datasheet: STAR136

BATCH NUMBER 154326

Description:	GOAT ANTI MOUSE IgG3
Specificity:	IgG3
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to mouse IgG3 were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.09% Sodium Azide

Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Mouse IgG3 paraproteins.
External Database Links	<p>UniProt: P03987 Related reagents</p> <p>Entrez Gene: 380795 AI324046 Related reagents</p>
RRID	AB_1102669
Specificity	Goat anti Mouse IgG3 antibody recognizes Mouse IgG3 and has been cross absorbed against mouse IgM, IgG1, IgG2a, IgG2b and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG3 antibody shows minimal cross-reactivity with human immunoglobulins.
References	<ol style="list-style-type: none"> Knipping, K. <i>et al.</i> (2011) A gastrointestinal rotavirus infection mouse model for immune modulation studies. Virology 8: 109. Donius LR <i>et al.</i> (2013) Optimal germinal center B cell activation and T-dependent antibody responses require expression of the mouse complement receptor Cr1. J Immunol. 191 (1): 434-47. Hwang, S.R. <i>et al.</i> (2015) Altered expression levels of neurodevelopmental proteins in fetal brains of BTBR T+tf/J mice with autism-like behavioral characteristics. J Toxicol Environ Health A. 78 (8): 516-23. Zhao, Z. <i>et al.</i> (2015) Multiple B-cell epitope vaccine induces a <i>Staphylococcus</i> enterotoxin B-specific IgG1 protective response against MRSA infection. Sci Rep. 5: 12371. Thema, N. <i>et al.</i> (2019) Identification and characterisation of conserved epitopes of <i>E. ruminantium</i> that activate Th1 CD4⁺ T cells: Towards the development of a multi-epitope vaccine. Mol Immunol. 107: 106-114.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10077 available at: https://www.bio-rad-antibodies.com/SDS/STAR13610077
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M369687:200529'

Printed on 25 Sep 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)