

Datasheet: STAR136P

Description:	GOAT ANTI MOUSE IgG3:HRP
Specificity:	IgG3
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA	▪			1/4000 - 1/8000
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Antisera to mouse IgG3 were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% Glycerol
Immunogen	Mouse IgG3 paraproteins.

External Database Links

UniProt:

[P03987](#) [Related reagents](#)

Entrez Gene:

[380795](#) [AI324046](#) [Related reagents](#)

RRID	AB_1102670
Specificity	Goat anti Mouse IgG3 antibody recognizes Mouse IgG3 and has been cross absorbed against mouse IgM, IgG1, IgG2a, IgG2b and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG3 antibody shows minimal cross-reactivity with human immunoglobulins.
References	<ol style="list-style-type: none"> Knipping, K. <i>et al.</i> (2011) A gastrointestinal rotavirus infection mouse model for immune modulation studies. Virology J. 8: 109. Donius LR <i>et al.</i> (2013) Optimal germinal center B cell activation and T-dependent antibody responses require expression of the mouse complement receptor Cr1. J Immunol. 191 (1): 434-47. Hwang, S.R. <i>et al.</i> (2015) Altered expression levels of neurodevelopmental proteins in fetal brains of BTBR T+tf/J mice with autism-like behavioral characteristics. J Toxicol Environ Health A. 78 (8): 516-23. Zhao, Z. <i>et al.</i> (2015) Multiple B-cell epitope vaccine induces a <i>Staphylococcus</i> enterotoxin B-specific IgG1 protective response against MRSA infection. Sci Rep. 5: 12371. Thema, N. <i>et al.</i> (2019) Identification and characterisation of conserved epitopes of <i>E. ruminantium</i> that activate Th1 CD4⁺ T cells: Towards the development of a multi-epitope vaccine. Mol Immunol. 107: 106-114.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10106 available at: 10106: https://www.bio-rad-antibodies.com/uploads/MSDS/10106.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M350478:190307'

Printed on 11 Aug 2020