

Datasheet: STAR138A

Description:	GOAT ANTI MOUSE IgM:Alk. Phos.(HUMAN ADSORBED)
Specificity:	IgM
Format:	Alk. Phos.
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA	▪			1/2000 - 1/4000
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG conjugated to alkaline phosphatase - liquid
Antiserum Preparation	Antisera to mouse IgM were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on IgM.
Buffer Solution	TRIS buffered saline, 1mM MgCl ₂
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃) 50% Glycerol
Immunogen	Mouse IgM paraproteins.

**External Database
Links**

UniProt:

[P01872](#) [Related reagents](#)
[P01873](#) [Related reagents](#)

Entrez Gene:

[16019](#) Ighm [Related reagents](#)
[16019](#) Ighm [Related reagents](#)

RRID AB_2076071

Specificity **Goat anti Mouse IgM antibody** recognizes mouse IgM. Goat anti Mouse IgM antibody has been cross absorbed against mouse IgA, IgG1, IgG2a, IgG2b and IgG3, pooled human sera and purified human paraproteins. Goat anti Mouse IgM antibody shows minimal cross reactivity with human immunoglobulins.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee Guaranteed until date of expiry. Please see product label.

**Health And Safety
Information** Material Safety Datasheet documentation #10322 available at:
10322: <https://www.bio-rad-antibodies.com/uploads/MSDS/10322.pdf>

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
--------------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M387611:210629'

Printed on 29 Aug 2021