

Datasheet: STAR141P

Description:	GOAT ANTI HUMAN IgA:HRP
Specificity:	IgA
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin	▪			
ELISA	▪			1/4000 - 1/8000
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Antisera to human IgA were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on pooled human IgA covalently linked to agarose.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative	0.01% Thiomersal
Stabilisers	50% Glycerol

Immunogen	Human IgA paraproteins.
------------------	-------------------------

External Database Links	UniProt:
--------------------------------	-----------------

[P01876](#) [Related reagents](#)

[P01877](#) [Related reagents](#)

Entrez Gene:

[3493](#) IGHA1 [Related reagents](#)

[3494](#) IGHA2 [Related reagents](#)

RRID AB_1102621

Specificity **Goat anti human IgA polyclonal antibody** recognizes the heavy chain of human IgA. This anti human IgA antibody has been cross adsorbed against human IgM and IgG but might cross-react with IgA from other species.

References 1. Whelan, S.F. *et al.* (2013) Distinct characteristics of antibody responses against factor VIII in healthy individuals and in different cohorts of hemophilia A patients. [Blood. 121: 1039-48.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10097 available at: 10097: <https://www.bio-rad-antibodies.com/uploads/MSDS/10097.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M381883:210512'

Printed on 21 Mar 2022
