

Datasheet: STAR144B

Description:	GOAT F(ab') ₂ ANTI HUMAN IgD HEAVY CHAIN:Biotin
Specificity:	IgD HEAVY CHAIN
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	F(ab') ₂ fragment of purified IgG conjugated to Biotin - liquid

Antiserum Preparation Antisera to human IgD were raised by repeated immunisation of goats with purified antigen. F(ab')₂ fragments were prepared by gel filtration chromatography of pepsin digested goat anti human IgD antibody.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
---------------------------------	--

Approx. Protein Concentrations	F(ab') ₂ concentration 0.5mg/ml
---------------------------------------	--

Immunogen Human IgD.

External Database

Links

UniProt:

[P01880](#) [Related reagents](#)

Entrez Gene:

[3495](#) IGHD [Related reagents](#)

RRID AB_10612766

Specificity **Goat F(ab')₂ anti Human IgD Heavy Chain antibody** recognizes the heavy chain of human IgD and has been cross absorbed against human IgM, IgA and IgG. The antibody might cross react with mouse IgD.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/STAR144B>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

IgG (H/L) (602001...) [Alk. Phos.](#)

IgG (H/L) (642001...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M428712:240301'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)