

Datasheet: STAR146

BATCH NUMBER 153874

Description:	GOAT F(ab') ₂ ANTI HUMAN IgM
Specificity:	IgM
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	F(ab') ₂ fragments of purified IgG - liquid
Antiserum Preparation	Antisera to human IgM were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography. F(ab') ₂ fragments were prepared by pepsin digestion.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	F(ab') ₂ concentration 0.5 mg/ml

Immunogen Human IgM.

External Database

Links

UniProt:

[P01871](#) [Related reagents](#)

Entrez Gene:

[3507](#) IGHM [Related reagents](#)

RRID

AB_1102675

Specificity

Goat F(ab')₂ anti Human IgM antibody recognizes the heavy chain of human IgM and has been cross absorbed against human IgA and IgG. Goat F(ab')₂ anti Human IgM antibody might cross react with IgM from other species.

References

1. Mpakou, V.E. *et al.* (2017) Quantitative and qualitative analysis of regulatory T cells in B cell chronic lymphocytic leukemia. [Leuk Res. 60: 74-81.](#)
 2. Lu, D.R. *et al.* (2014) Identifying functional anti-*Staphylococcus aureus* antibodies by sequencing antibody repertoires of patient plasmablasts. [Clin Immunol. 152 \(1-2\): 77-89.](#)
 3. Tan, Y.C. *et al.* (2014) High-throughput sequencing of natively paired antibody chains provides evidence for original antigenic sin shaping the antibody response to influenza vaccination. [Clin Immunol. 151 \(1\): 55-65.](#)
-

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information

Material Safety Datasheet documentation #10123 available at:
<https://www.bio-rad-antibodies.com/SDS/STAR146>
10123

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350484:190307'

Printed on 28 Aug 2024

