

## Datasheet: STAR146F

**BATCH NUMBER 155249**

<b>Description:</b>	GOAT F(ab') <sub>2</sub> ANTI HUMAN IgM:FITC
<b>Specificity:</b>	IgM
<b>Format:</b>	FITC
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
Immunofluorescence	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

#### Target Species

Human

#### Product Form

F(ab')<sub>2</sub> fragments of purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525

#### Antiserum Preparation

Antisera to human IgM were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography. F(ab')<sub>2</sub> fragments were prepared by pepsin digestion.

#### Buffer Solution

Phosphate buffered saline

#### Preservative Stabilisers

0.1% Sodium Azide (NaN<sub>3</sub>)

<b>Approx. Protein Concentrations</b>	F(ab') <sub>2</sub> concentration 0.5 mg/ml
<b>Immunogen</b>	Human IgM.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P01871</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">3507</a>    IGHM    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_1102673
<b>Specificity</b>	<b>Goat F(ab')<sub>2</sub> anti Human IgM antibody</b> recognizes the heavy chain of human IgM and has been cross absorbed against human IgA and IgG. Goat F(ab') <sub>2</sub> anti Human IgM antibody might cross react with IgM from other species.
<b>Flow Cytometry</b>	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells in 100ul.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Mpakou, V.E. <i>et al.</i> (2017) Quantitative and qualitative analysis of regulatory T cells in B cell chronic lymphocytic leukemia. <a href="#">Leuk Res. 60: 74-81.</a></li> <li>2. Lu, D.R. <i>et al.</i> (2014) Identifying functional anti-<i>Staphylococcus aureus</i> antibodies by sequencing antibody repertoires of patient plasmablasts. <a href="#">Clin Immunol. 152 (1-2): 77-89.</a></li> <li>3. Tan, Y.C. <i>et al.</i> (2014) High-throughput sequencing of natively paired antibody chains provides evidence for original antigenic sin shaping the antibody response to influenza vaccination. <a href="#">Clin Immunol. 151 (1): 55-65.</a></li> </ol>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/STAR146F10040">https://www.bio-rad-antibodies.com/SDS/STAR146F10040</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M358117:190911'

**Printed on 01 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)