

## Datasheet: STAR147A

**BATCH NUMBER 170788**

<b>Description:</b>	GOAT ANTI HUMAN IgE:Alk. Phos.
<b>Specificity:</b>	IgE
<b>Format:</b>	Alk. Phos.
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/1000 - 1/10000
Immunoprecipitation			▪	
Western Blotting	▪			1/1000 - 1/5000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG conjugated to Alkaline Phosphatase - liquid

**Antiserum Preparation** Antiserum to human IgE was raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

**Buffer Solution** 50mM HEPES, 0.1M NaCl, 1mM MgCl<sub>2</sub>, 0.1mM ZnCl<sub>2</sub>

**Preservative Stabilisers**  
0.09% Sodium Azide  
0.2% Bovine Serum Albumin

**Approx. Protein Concentrations** IgG concentration 0.5 mg/ml.

<b>Immunogen</b>	Native human IgE.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P01854</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">3497</a>    IGHE    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_1221598
<b>Specificity</b>	<p><b>Goat anti Human IgE polyclonal antibody</b> recognizes human IgE. Minimal cross-reactivity (&lt;0.1%) is expected with other human immunoglobulins and light chains.</p> <p>Goat anti Human IgE in a biotinylated format has been successfully employed as a detection reagent in the development of a sensitive sandwich ELISA for the measurement of human IgE.</p>
<b>References</b>	<ol style="list-style-type: none"> <li>1. Khan, F.M. <i>et al.</i> (2012) Basophil activation test compared to skin prick test and fluorescence enzyme immunoassay for aeroallergen-specific Immunoglobulin-E. <a href="#">Allergy Asthma Clin Immunol. 8 (1): 1.</a></li> <li>2. Pinot de Moira, A. <i>et al.</i> (2013) Effects of treatment on IgE responses against parasite allergen-like proteins and immunity to reinfection in childhood schistosome and hookworm coinfections. <a href="#">Infect Immun. 81 (1): 23-32.</a></li> <li>3. Zheng, W. <i>et al.</i> (2016) Upregulated expression of substance P in basophils of the patients with chronic spontaneous urticaria: induction of histamine release and basophil accumulation by substance P. <a href="#">Cell Biol Toxicol. 32 (3): 217-28.</a></li> <li>4. Jin Y <i>et al.</i> (2015) Allergic response to squid (<i>Todarodes pacificus</i>) tropomyosin Tod p1 structure modifications induced by high hydrostatic pressure. <a href="#">Food Chem Toxicol. 76: 86-93.</a></li> <li>5. Wei, J.F. <i>et al.</i> (2009) Induction of mast cell accumulation, histamine release and skin edema by N49 phospholipase A2. <a href="#">BMC Immunol. 10: 21.</a></li> <li>6. Tsai, C.T. <i>et al.</i> (2018) Isotype-specific agglutination-PCR (ISAP): A sensitive and multiplex method for measuring allergen-specific IgE. <a href="#">J Allergy Clin Immunol. 141 (5): 1901-4.e15.</a></li> </ol>
<b>Storage</b>	<p>Store at +4°C.  DO NOT FREEZE</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch.
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10089 available at:  <a href="https://www.bio-rad-antibodies.com/SDS/STAR147A">https://www.bio-rad-antibodies.com/SDS/STAR147A</a>  10089</p>
<b>Regulatory</b>	For research purposes only.

**North & South** Tel: +1 800 265 7376

**Worldwide** Tel: +44 (0)1865 852 700

**Europe** Tel: +49 (0) 89 8090 95 21

**America** Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M428720:240301'

**Printed on 23 May 2025**

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)