

## Datasheet: STAR165

**BATCH NUMBER 172793**

<b>Description:</b>	GOAT ANTI HUMAN LAMBDA LIGHT CHAIN (MOUSE ADSORBED)
<b>Specificity:</b>	LAMBDA LIGHT CHAIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting	▪			
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**Target Species** Human

**Product Form** Purified IgG - liquid

**Antiserum Preparation** Antisera to human lambda chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

**Buffer Solution** Borate buffered saline

**Preservative Stabilisers** <0.1% Sodium Azide (NaN<sub>3</sub>)

**Approx. Protein Concentrations** IgG concentration 1.0mg/ml

<b>Immunogen</b>	Human myeloma proteins with lambda light chains.
<b>RRID</b>	AB_10601744
<b>Specificity</b>	<b>Goat anti Human lambda light chain antibody</b> recognises human lambda light chains. Goat anti Human lambda light chain antibody has been cross adsorbed against pooled human myeloma proteins with kappa light chains and mouse immunoglobulins.
<b>Storage</b>	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.  Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10077 available at: <a href="https://www.bio-rad-antibodies.com/SDS/STAR165">https://www.bio-rad-antibodies.com/SDS/STAR165</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M428732:240301'

Printed on 29 Jan 2026