

## Datasheet: STAR16D488GA

**BATCH NUMBER 0911**

<b>Description:</b>	RABBIT F(ab') <sub>2</sub> ANTI RAT IgG:DyLight®488
<b>Specificity:</b>	IgG
<b>Format:</b>	DyLight®488
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/5 - 1/400
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls

<b>Target Species</b>	Rat		
<b>Product Form</b>	Purified IgG conjugated to DyLight®488 - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	DyLight®488	493	518

**Antiserum Preparation** Antisera to rat IgG were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography

**Buffer Solution** Phosphate buffered saline

**Preservative Stabilisers** 0.09% Sodium Azide

**Approx. Protein Concentrations** IgG concentration 1.0 mg/ml

**Immunogen** Purified polyclonal rat IgG.

---

**External Database**

**Links**

**UniProt:**

[P20759](#) [Related reagents](#)  
[P20762](#) [Related reagents](#)  
[P20761](#) [Related reagents](#)  
[P20760](#) [Related reagents](#)

**Entrez Gene:**

[299354](#) Ighg [Related reagents](#)  
[362795](#) LOC362795 [Related reagents](#)  
[679045](#) LOC679045 [Related reagents](#)

---

**RRID** AB\_10846819

---

**Specificity** **Rabbit F(ab')<sub>2</sub> anti Rat IgG antibody** recognizes rat IgG, detecting all subclasses in immunodiffusion assays.

---

**Flow Cytometry** Use 50ul of the suggested working dilution to label 1x10<sup>6</sup> cells in 100ul

---

**Storage** Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/STAR16D488GA>  
10040

---

**Regulatory** For research purposes only

---

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M369726:200529'

**Printed on 18 Jan 2024**

---

