

## Datasheet: STAR208P

**BATCH NUMBER 166194**

<b>Description:</b>	GOAT ANTI RABBIT IgG (H/L):HRP
<b>Specificity:</b>	IgG (H/L)
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	2 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/5000 - 1/50,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit
<b>Product Form</b>	Purified IgG - lyophilized
<b>Reconstitution</b>	Reconstitute with 2.0 ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Antiserum Preparation</b>	Antisera to rabbit IgG (H/L) were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography using antigens coupled to agarose beads.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	15 mg/ml Bovine Serum Albumin
<b>Approx. Protein Concentrations</b>	IgG concentration 0.8 mg/ml after reconstitution

<b>Specificity</b>	<b>Goat anti Rabbit IgG (H/L) polyclonal antibody</b> reacts with the heavy chains of rabbit IgG, and the light chains common to the majority of rabbit immunoglobulins. No antibody reactivity was detected against non-immunoglobulin rabbit serum proteins. This antibody may cross-react with immunoglobulins from other species.
<b>Storage</b>	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10363 available at: <a href="https://www.bio-rad-antibodies.com/SDS/STAR208P10363">https://www.bio-rad-antibodies.com/SDS/STAR208P10363</a>
<b>Regulatory</b>	For research purposes only

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M421326:230706'

**Printed on 01 Mar 2024**