

## Datasheet: STAR3030P

**BATCH NUMBER 157323**

<b>Description:</b>	GOAT ANTI RAT IgG:HRP
<b>Specificity:</b>	IgG
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
ELISA	▪			1/4000 - 1/8000
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rat
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on pooled rat IgG linked to agarose
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	50% Glycerol
<b>Immunogen</b>	Normal rat IgG

#### Specificity

**Goat anti Rat IgG antibody** is specific for the gamma chain of rat IgG. It may also react with IgG from other species.

Goat anti Rat IgG antibody has been cross-adsorbed against rat IgM to minimize cross-reactivity.

---

<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10106 available at: <a href="https://www.bio-rad-antibodies.com/SDS/STAR3030P">https://www.bio-rad-antibodies.com/SDS/STAR3030P</a>
<b>Regulatory</b>	For research purposes only

---

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M353996:190620'

**Printed on 03 Sep 2025**

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)