

## Datasheet: STAR35A

**BATCH NUMBER 151294**

<b>Description:</b>	SHEEP ANTI RABBIT IgG:RPE
<b>Specificity:</b>	IgG
<b>Format:</b>	RPE
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/5 - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit		
<b>Product Form</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
<b>Reconstitution</b>	Reconstitute with 1 ml distilled water		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578

**Antiserum Preparation** Antisera to rabbit IgG were raised by repeated immunisation of sheep with highly purified antigen. Purified IgG prepared by affinity chromatography

**Buffer Solution** Phosphate buffered saline

**Preservative** 0.09% Sodium Azide  
**Stabilisers** 1% Bovine Serum Albumin  
 5% Sucrose

**Immunogen** Rabbit IgG.

---

**External Database Links**

**UniProt:**

[P01870](#)    [Related reagents](#)

**Entrez Gene:**

[100009097](#)    LOC100009097    [Related reagents](#)

---

**RRID** AB\_321858

---

**Specificity** **Sheep anti Rabbit IgG antibody** recognizes rabbit IgG.

---

**Flow Cytometry** Use 50ul of the suggested working dilution to label 10<sup>6</sup> cells in 100ul.

---

**References**

1. Båge, T. *et al.* (2010) Signal pathways JNK and NF-kappaB, identified by global gene expression profiling, are involved in regulation of TNFalpha-induced mPGES-1 and COX-2 expression in gingival fibroblasts. [BMC Genomics. 11: 241.](#)
  2. Senthilkumaran, C. *et al.* (2015) Localization of annexins A1 and A2 in the respiratory tract of healthy calves and those experimentally infected with *Mannheimia haemolytica*. [Vet Res. 46: 6.](#)
  3. Williams, K.M. *et al.* (2022) Prospective Phase II Trial of Montelukast to Treat Bronchiolitis Obliterans Syndrome after Hematopoietic Cell Transplantation and Investigation into Bronchiolitis Obliterans Syndrome Pathogenesis. [Transplant Cell Ther. S2666-6367\(22\)00049-5.](#)
- 

**Storage**

Prior to reconstitution store at +4°C. Following reconstitution store at +4°C.

DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #20487 available at: <https://www.bio-rad-antibodies.com/SDS/STAR35A>  
20487

---

**Regulatory**

For research purposes only

---

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M375710:210104'

**Printed on 19 Jan 2024**

