

Datasheet: STAR71D650

Description:	GOAT ANTI RAT IgG:DyLight®650 (MOUSE ADSORBED)
Specificity:	IgG (MOUSE ADSORBED)
Format:	DyLight®650
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/100 - 1/500
Immunofluorescence	▪			1/100 - 1/500

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat		
Product Form	Purified IgG conjugated to DyLight®650 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Dylight®650	654	673

Antiserum Preparation Antisera to rat IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG prepared by affinity chromatography

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Rat IgG.

**External Database
Links**

UniProt:

[P20759](#) [Related reagents](#)
[P20762](#) [Related reagents](#)
[P20761](#) [Related reagents](#)
[P20760](#) [Related reagents](#)

Entrez Gene:

[299354](#) Ighg [Related reagents](#)
[362795](#) LOC362795 [Related reagents](#)
[679045](#) LOC679045 [Related reagents](#)

Specificity **Goat anti Rat IgG (Mouse Adsorbed) antibody** recognizes rat IgG. Cross-reactivity with mouse IgG has been minimised by adsorption.

Flow Cytometry Use 50ul of the suggested working dilution to label 1x10⁶ cells in 100ul.

References

1. Yang, X. *et al.* (2010) The role of the JAK2-STAT3 pathway in pro-inflammatory responses of EMF-stimulated N9 microglial cells. [J Neuroinflammation. 7: 54.](#)
2. Tamayo, J. *et al.* (2001) Chemical sensors and biosensors in liquid environment based on microcantilevers with amplified quality factor. [Ultramicroscopy. 86: 167-73.](#)
3. Pérez-Bosque A *et al.* (2004) Dietary plasma protein affects the immune response of weaned rats challenged with *S. aureus* Superantigen B. [J Nutr. 134: 2667-72.](#)
4. Balan, P. *et al.* (2010) Immunomodulatory effects of ovine serum immunoglobulin in the growing rat. [Animal. 4: 1702-8.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

**Health And Safety
Information** Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
--------------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M390357:210908'

