

Datasheet: STAR71D800GA

Description:	GOAT ANTI RAT IgG:DyLight®800 (MOUSE ADSORBED)
Specificity:	IgG (MOUSE ADSORBED)
Format:	DyLight®800
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

RRID AB_10852249

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Western Blotting	▪			1/10000 - 1/50000
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls

Target Species Rat

Product Form Purified IgG conjugated to DyLight®800 - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	DyLight®800	777	794

Antiserum Preparation Antisera to rat IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG prepared by affinity chromatography

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen Rat IgG.

External Database Links
UniProt:
[P20759](#) [Related reagents](#)

[P20762](#) [Related reagents](#)

[P20761](#) [Related reagents](#)

[P20760](#) [Related reagents](#)

Entrez Gene:

[299354](#) Ighg [Related reagents](#)

[362795](#) LOC362795 [Related reagents](#)

[679045](#) LOC679045 [Related reagents](#)

Specificity **Goat anti Rat IgG (Mouse Adsorbed) antibody** recognizes rat IgG. Cross-reactivity with mouse IgG has been minimised by adsorption.

Flow Cytometry Use 50ul of the suggested working dilution to label 1x10⁶ cells in 100ul

References

1. Yang, X. *et al.* (2010) The role of the JAK2-STAT3 pathway in pro-inflammatory responses of EMF-stimulated N9 microglial cells. [J Neuroinflammation. 7: 54.](#)
2. Tamayo, J. *et al.* (2001) Chemical sensors and biosensors in liquid environment based on microcantilevers with amplified quality factor. [Ultramicroscopy. 86: 167-73.](#)
3. Pérez-Bosque A *et al.* (2004) Dietary plasma protein affects the immune response of weaned rats challenged with *S. aureus* Superantigen B. [J Nutr. 134: 2667-72.](#)
4. Balan, P. *et al.* (2010) Immunomodulatory effects of ovine serum immunoglobulin in the growing rat. [Animal. 4: 1702-8.](#)

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use

Guarantee 18 months from date of despatch

Acknowledgements DyLight[®] is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M353031:190409'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)