

Datasheet: STAR74D800GA

Description:	GOAT ANTI MOUSE IgG:DyLight®800 (RAT ADSORBED)
Specificity:	IgG (RAT ADSORBED)
Format:	DyLight®800
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

RRID AB_10846934

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Western Blotting	▪			1/5000 - 1/20000
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls

Target Species Mouse

Product Form Purified IgG conjugated to DyLight®800 - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Dylight®800	777	794

Antiserum Preparation Antisera to mouse IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG prepared by affinity chromatography

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen Mouse IgG.

External Database Links **UniProt:**

[P01869](#) [Related reagents](#)
[P03987](#) [Related reagents](#)
[P01867](#) [Related reagents](#)
[P01864](#) [Related reagents](#)
[P01868](#) [Related reagents](#)
[P01863](#) [Related reagents](#)

Entrez Gene:

[16017](#) Ighg1 [Related reagents](#)
[380795](#) AI324046 [Related reagents](#)
[16016](#) Ighg2b [Related reagents](#)
[16017](#) Ighg1 [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)

Synonyms	Igh-4
-----------------	-------

Specificity	Goat anti Mouse IgG (Rat Adsorbed) antibody recognizes mouse IgG. Cross-reactivity with rat IgG has been minimised by adsorption.
--------------------	--

Flow Cytometry	Use 50ul of the suggested working dilution to label 1×10^6 cells in 100ul
-----------------------	--

References	1. Jarrett, B.R. <i>et al.</i> (2010) <i>In vivo</i> mapping of vascular inflammation using multimodal imaging. PLoS One. 5(10). pii: e13254.
-------------------	---

Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use
----------------	---

Guarantee	18 months from date of despatch
------------------	---------------------------------

Acknowledgements	DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries
-------------------------	---

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M353038:190409'

Printed on 20 May 2019