

## Datasheet: STAR79P

<b>Description:</b>	GOAT ANTI HAMSTER IgG:HRP
<b>Specificity:</b>	IgG
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.8 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/1000 - 1/5000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Hamster
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

**Antiserum Preparation** Antisera to hamster IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography.

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal 1% Bovine Serum Albumin

<b>Approx. Protein Concentrations</b>	IgG concentration 0.8 mg/ml
---------------------------------------	-----------------------------

<b>Immunogen</b>	Hamster IgG.
------------------	--------------

<b>RRID</b>	AB_321907
-------------	-----------

<b>Specificity</b>	<b>Goat anti Hamster IgG antibody</b> recognizes Golden Syrian and Armenian hamster IgG (H+L) and has been adsorbed against both mouse and rat immunoglobulins to minimise cross-reactivity.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Moore, G.T. <i>et al.</i> (2008) Glycosylation changes in hFUT1 transgenic mice increase TCR signaling and apoptosis resulting in thymocyte maturation arrest. <a href="#">Mol Immunol. 45: 2401-10.</a></li> <li>2. Samant, M. <i>et al.</i> (2009) Immunization with the DNA-encoding N-terminal domain of proteophosphoglycan of <i>Leishmania donovani</i> generates Th1-type immunoprotective response against experimental visceral leishmaniasis. <a href="#">J Immunol. 183: 470-9.</a></li> <li>3. Bacelo, K.L. <i>et al.</i> (2014) Xanthan gum as an adjuvant in a subunit vaccine preparation against leptospirosis. <a href="#">Biomed Res Int. 2014: 636491.</a></li> <li>4. Verma R <i>et al.</i> (2015) Cross reactive molecules of human lymphatic filaria <i>Brugia malayi</i> inhibit <i>Leishmania donovani</i> infection in hamsters. <a href="#">Acta Trop. 152: 103-11.</a></li> <li>5. Forster, K.M. <i>et al.</i> (2015) DNA prime-protein boost based vaccination with a conserved region of leptospiral immunoglobulin-like A and B proteins enhances protection against leptospirosis. <a href="#">Mem Inst Oswaldo Cruz. 110 (8): 989-95.</a></li> <li>6. Wiśniewski, M. <i>et al.</i> (2016) Hamsters vaccinated with Ace-mep-7 DNA vaccine produced protective immunity against <i>Ancylostoma ceylanicum</i> infection. <a href="#">Exp Parasitol. 163: 1-7.</a></li> <li>7. Verma, R. <i>et al.</i> (2018) <i>Leishmania donovani</i>. molecules recognized by sera of filaria infected host facilitate filarial infection. <a href="#">Parasitol Res. 117 (9): 2901-12.</a></li> <li>8. Spitzova, T. <i>et al.</i> (2020) Interactions between host biogenic amines and sand fly salivary yellow-related proteins. <a href="#">Parasit Vectors. 13 (1): 237.</a></li> <li>9. Fernández, L. <i>et al.</i> (2021) Protective Efficacy in a Hamster Model of a Multivalent Vaccine for Human Visceral Leishmaniasis (MuLeVaClin) Consisting of the KMP11, LEISH-F3+, and LJL143 Antigens in Virosomes, Plus GLA-SE Adjuvant. <a href="#">Microorganisms. 9 (11): 2253.</a></li> </ol>
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10095 available at: <a href="https://www.bio-rad-antibodies.com/SDS/STAR79P">https://www.bio-rad-antibodies.com/SDS/STAR79P</a> 10095
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)  
[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)  
[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)  
[TMB CORE \(BUF056A\)](#)  
[TMB CORE+ \(BUF062A\)](#)  
[TMB SIGNAL+ \(BUF054A\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M428764:240301'

**Printed on 20 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)