

Datasheet: STAR87A

BATCH NUMBER 150602

| | |
|----------------------|------------------------------------|
| Description: | GOAT ANTI MOUSE IgG/A/M:Alk. Phos. |
| Specificity: | IgG IgA IgM |
| Format: | Alk. Phos. |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Immunohistology - Frozen | ▪ | | | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/3000 |
| Western Blotting | ▪ | | | |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Product Form

Purified IgG conjugated to alkaline phosphatase - liquid

Antiserum Preparation

Antisera to mouse immunoglobulins were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

0.05M TRIS Chloride
0.15M NaCl
0.001M MgCl₂
0.0001M ZnCl₂
50% (v/v) Glycerol; pH8.0

Preservative

0.1% Sodium Azide (NaN₃)

Stabilisers

1% Bovine Serum Albumin

| | |
|---------------------------------------|--|
| Approx. Protein Concentrations | IgG concentration 1.0 mg/ml |
| Immunogen | Mouse IgG purified from normal mouse serum, mouse IgA and IgM purified from ascitic fluid. |
| External Database Links | <p>UniProt:</p> <p>P01867 Related reagents</p> <p>P01865 Related reagents</p> <p>P01864 Related reagents</p> <p>P01863 Related reagents</p> <p>P01869 Related reagents</p> <p>P01868 Related reagents</p> <p>P03987 Related reagents</p> <p>P01872 Related reagents</p> <p>P01873 Related reagents</p> <p>P01878 Related reagents</p> <p>Entrez Gene:</p> <p>16016 Ighg2b Related reagents</p> <p>380793 Igh-1a Related reagents</p> <p>16017 Ighg1 Related reagents</p> <p>16017 Ighg1 Related reagents</p> <p>380793 Igh-1a Related reagents</p> <p>380793 Igh-1a Related reagents</p> <p>16019 Ighm Related reagents</p> <p>16019 Ighm Related reagents</p> <p>16061 Igh-VJ558 Related reagents</p> <p>380795 AI324046 Related reagents</p> |
| Synonyms | Igh-4 |
| RRID | AB_321852 |
| Specificity | Goat anti Mouse IgG/A/M (Human Adsorbed) antibody recognizes all subclasses of murine immunoglobulin and has been adsorbed against human, bovine and hamster immunoglobulins. |
| References | <ol style="list-style-type: none"> 1. Bergmeier, L.A. <i>et al.</i> (2005) Mucosal alloimmunization elicits T-cell proliferation, CC chemokines, CCR5 antibodies and inhibition of simian immunodeficiency virus infectivity. J Gen Virol. 86: 2231-8. 2. Peters, B. <i>et al.</i> (2004) Effect of heterosexual intercourse on mucosal alloimmunisation and resistance to HIV-1 infection. Lancet. 363: 518-24. 3. Bartlomiejczyk, M.A. <i>et al.</i> (2014) Interaction of lectin pathway of complement-activating pattern recognition molecules with mycobacteria. Clin Exp Immunol. 178 (2): 310-9. |

4. Kubelkova, K. *et al.* (2021) Early infection-induced natural antibody response. [Sci Rep. 11 \(1\): 1541.](#)
5. Grandoni, F. *et al.* (2017) Characterization of leukocyte subsets in buffalo (*Bubalus bubalis.*) with cross-reactive monoclonal antibodies specific for bovine MHC class I and class II molecules and leukocyte differentiation molecules. [Dev Comp Immunol. 74: 101-109.](#)

Storage

Store at +4°C.

DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10063 available at: <https://www.bio-rad-antibodies.com/SDS/STAR87A>
10063

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M415437:221222'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)