

Datasheet: STAR87F

Description:	GOAT ANTI MOUSE IgG/A/M:FITC (HUMAN ADSORBED)
Specificity:	IgG IgA IgM
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1.8 mg

Product Details

RRID AB_321853

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/100 - 1/500
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Mouse

Product Form Purified IgG - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525

Antiserum Preparation Antisera to mouse Ig were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 1.8 mg/ml

Immunogen Mouse IgG purified from normal mouse serum, mouse IgA and IgM purified from ascitic fluid.

External Database Links **UniProt:**

[P01867](#) [Related reagents](#)
[P01865](#) [Related reagents](#)
[P01864](#) [Related reagents](#)
[P01863](#) [Related reagents](#)
[P01869](#) [Related reagents](#)
[P01868](#) [Related reagents](#)
[P03987](#) [Related reagents](#)
[P01872](#) [Related reagents](#)
[P01873](#) [Related reagents](#)
[P01878](#) [Related reagents](#)

Entrez Gene:

[16016](#) Ighg2b [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[16017](#) Ighg1 [Related reagents](#)
[16017](#) Ighg1 [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[16019](#) Ighm [Related reagents](#)
[16019](#) Ighm [Related reagents](#)
[16061](#) Igh-VJ558 [Related reagents](#)
[380795](#) AI324046 [Related reagents](#)

Synonyms

Igh-4

Specificity

Goat anti Mouse IgG/A/M (Human Adsorbed) antibody recognizes all subclasses of murine immunoglobulin and has been adsorbed against human, bovine and hamster immunoglobulins.

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Bergmeier, L.A. *et al.* (2005) Mucosal alloimmunization elicits T-cell proliferation, CC chemokines, CCR5 antibodies and inhibition of simian immunodeficiency virus infectivity. [J Gen Virol. 86: 2231-8.](#)
2. Peters, B. *et al.* (2004) Effect of heterosexual intercourse on mucosal alloimmunisation and resistance to HIV-1 infection. [Lancet. 363: 518-24.](#)
3. Bartlomiejczyk, M.A. *et al.* (2014) Interaction of lectin pathway of complement-activating pattern recognition molecules with mycobacteria. [Clin Exp Immunol. 178 \(2\): 310-9.](#)

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

18 months from date of despatch.

**Health And Safety
Information**

Material Safety Datasheet documentation #10041 available at:
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M353047:190409'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)