

Datasheet: STAR8D405GA

BATCH NUMBER 0911

Description:	RABBIT F(ab') ₂ ANTI MOUSE IgG:DyLight®405
Specificity:	IgG
Format:	DyLight®405
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/20 - 1/200
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls

Target Species	Mouse								
Species Cross Reactivity	Reacts with: Rat N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.								
Product Form	Purified IgG conjugated to DyLight®405 - liquid								
Max Ex/Em	<table><tr><th>Fluorophore</th><th>Excitation Max (nm)</th><th>Emission Max (nm)</th></tr><tr><td>Dylight®405</td><td>400</td><td>420</td></tr></table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	Dylight®405	400	420		
Fluorophore	Excitation Max (nm)	Emission Max (nm)							
Dylight®405	400	420							
Preparation	Purified IgG was prepared from whole serum by affinity chromatography. F(ab')2 fragments were prepared by pepsin digestion of the IgG followed by a gel filtration step to remove any remaining intact IgG or Fc fragments								
Buffer Solution	Phosphate buffered saline								

Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Whole mouse IgG.
External Database Links	UniProt: P01869 Related reagents P01865 Related reagents P03987 Related reagents P01867 Related reagents P01868 Related reagents P01864 Related reagents P01863 Related reagents Entrez Gene: 16017 Ighg1 Related reagents 380793 Igh-1a Related reagents 16016 Ighg2b Related reagents 16017 Ighg1 Related reagents 380793 Igh-1a Related reagents 380795 AI324046 Related reagents 380793 Igh-1a Related reagents
Synonyms	Igh-4
RRID	AB_10843452
Specificity	Rabbit F(ab')₂ anti Mouse IgG antibody recognizes all subclasses of mouse IgG. Cross reactivity with rat IgG is expected. Cross reactivity with human serum proteins has been minimised by solid phase adsorption. STAR8B is suitable as a bridge in a PAP / APAAP complex. Peroxidase anti-Peroxidase (PAP) complexes are formed by antibodies raised towards Horse Radish Peroxidase (HRP), binding to HRP molecules. Such complexes are used as signal amplification reagents in immunochemical staining techniques. It must be emphasised that the primary antibody must be raised in the same species as the antibodies used to form the PAP complex, so that a secondary "linker" antibody can be used to "bridge" the primary antibody to the PAP complex.
Flow Cytometry	Use 50ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul
References	1. Seehafer, S.S. <i>et al.</i> (2011) Immunosuppression alters disease severity in juvenile Batten disease mice. J Neuroimmunol. 230 (1-2): 169-72.

2. Motallebzadeh, R. *et al.* (2012) Blocking lymphotoxin signaling abrogates the development of ectopic lymphoid tissue within cardiac allografts and inhibits effector antibody responses. [FASEB J. 26 \(1\): 51-62.](#)
3. Niven, J. *et al.* (2015) S100B Up-Regulates Macrophage Production of IL1 β and CCL22 and Influences Severity of Retinal Inflammation. [PLoS One. 10 \(7\): e0132688.](#)

Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use
Guarantee	12 months from date of despatch
Acknowledgements	DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/STAR8D405GA 10040
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
-----------------------	--	-----------	--	--------	--

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M369802:200529'

Printed on 18 Jan 2024