

Datasheet: TZA003

BATCH NUMBER 157944

Description:	HuCAL Fab-F-Spy2-H NEGATIVE CONTROL ANTIBODY
Specificity:	HuCAL Fab-F-Spy2-H NEGATIVE CONTROL
Format:	Purified
Product Type:	Negative/Isotype Control
Clone:	AbD18705ad
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Negative Control
Product Form	A monovalent human recombinant Fab (kappa light chain) selected from the HuCAL® phage display library. Expressed in a proprietary <i>E. coli</i> strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as liquid.
Preparation	Metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% ProClin 950
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Green fluorescent protein

Specificity **HuCAL Fab-F-Spy2-H negative control, clone AbD18705ad** is a recombinant Fab antibody fragment in the format Fab-F-Spy2-H with specificity for green fluorescent protein (GFP). It has no known reactivity with mammalian proteins or other antigens. It has been tested by ELISA for specificity to GFP, but has not been tested in other applications. It is suggested as a negative control reagent when using other HuCAL antibodies of the same Fab format.

It is recommended that this reagent is used at the same concentration as the test reagent.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements ProClin is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow
Sold under license of patent and patent applications, including but not limited to the international PCT application PCT/GB2018/051065, the European patent EP2534484B1, UK patent GB1706430.4, the US patents US9,547,003, US10,247,727, US patent application US16/274908 and corresponding patents
HuCAL® is a registered trademark of MorphoSys AG

His-tag is a registered trademark of EMD Biosciences.

Health And Safety Information Material Safety Datasheet documentation #20479 available at: <https://www.bio-rad-antibodies.com/SDS/TZA003>
20479

Licensed Use For *in vitro* research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at [SpyTag and SpyCatcher Products](#), and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[SpyCatcher2 \(TZC001\)](#)

[SpyCatcher2-CYS \(TZC001CYS\)](#)

[BiSpyCatcher2 \(TZC002\)](#)

[BiSpyCatcher2-CYS \(TZC002CYS\)](#)

[BiSpyCatcher2-CYS3 \(TZC002CYS3\)](#)

[HUMAN IgG1-FcSpyCatcher3 \(TZC009\)](#)
[MOUSE IgG2a-FcSpyCatcher3 \(TZC012\)](#)
[RABBIT IgG-FcSpyCatcher3 \(TZC013\)](#)
[HUMAN IgA-FcSpyCatcher3 \(TZC020\)](#)
[SpyCatcher3 \(TZC025\)](#)
[SpyCatcher3-CYS \(TZC025CYS\)](#)
[BiSpyCatcher2:HRP \(TZC002P\)](#)
[HUMAN IgG2-FcSpyCatcher3 \(TZC016\)](#)
[HUMAN IgG3-FcSpyCatcher3 \(TZC017\)](#)
[HUMAN IgG4-FcSpyCatcher3 \(TZC018\)](#)
[HUMAN IgG4-Pro-FcSpyCatcher3 \(TZC019\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M389940:210816'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)