

Datasheet: TZA009P

| | |
|----------------------|---|
| Description: | HUMAN ANTI INFLIXIMAB (DRUG/TARGET COMPLEX):HRP |
| Specificity: | INFLIXIMAB DRUG/TARGET COMPLEX |
| Other names: | Remicade |
| Format: | HRP |
| Product Type: | Monoclonal Antibody |
| Clone: | AbD18391 |
| Isotype: | Fab antibody |
| Quantity: | 0.1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A bivalent human recombinant Fab (lambda light chain) generated by coupling of the monovalent SpyTagged Fab [TZA009](#) with a HRP conjugated BiSpyCatcher version 2 ([TZC002P](#)). The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

| | |
|-------------------------|--|
| Immunogen | Infliximab/TNF α complex |
| Specificity | <p>Human Anti-Infliximab (Drug/Target Complex) Antibody, clone AbD18391 specifically recognizes the antibody drug infliximab when in complex with its target, TNFα (cytokine “tumour necrosis factor-alpha”). At defined experimental levels, the antibody does not recognize free infliximab nor unbound recombinant human TNFα. However, as levels greatly increase, detection is noted. It can be used to measure levels of infliximab in patient samples.</p> <p>This TrailBlazer Antibody is available in two different formats, monovalent Fab unconjugated (TZA009) and bivalent Fab conjugated to HRP (TZA009P). The monovalent recombinant Fab antibody has a SpyTag incorporated, which is used to create the bivalent format by coupling two SpyTagged Fab antibodies to an HRP conjugated BiSpyCatcher2 (TZC002P). Either format can be used to develop a pharmacokinetic (PK) antigen capture assay to measure free infliximab captured via immobilized TNFα.</p> <p>The antibody is tagged with a SpyTag2 at the C-terminus of the Fab antibody heavy chain, enabling the user to couple it to any of the SpyCatchers available in our catalog for conversion to alternative formats in less than one hour.</p> <p>View a summary of all SpyCatcher products</p> <p>Infliximab (Remicade) is a chimeric IgG1/kappa antibody. It has been approved for the treatment of psoriasis, Crohn’s Disease, ankylosing spondylitis, psoriatic arthritis, rheumatoid arthritis and ulcerative colitis. Infliximab acts by blocking the binding of TNFα to TNF receptors. Infliximab also causes programmed cell death of TNFα-expressing activated T lymphocytes.</p> <p>View a summary of all anti-infliximab antibodies</p> |
| Affinity | The monovalent intrinsic affinity of AbD18391 was measured as $K_D = 16$ nM by real time, label free molecular interaction analysis on immobilized biotinylated infliximab complexed with human TNF α . |
| ELISA | <p>Clone AbD18391pap can be used in a direct ELISA to detect the infliximab/TNFα complex or as the detection antibody for infliximab captured on immobilized TNFα in a PK antigen capture ELISA.</p> <p>Protocol: PK antigen capture ELISA to measure bound drug exclusively.</p> |
| Storage | <p>This product is shipped frozen.</p> <p>When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.</p> |
| Guarantee | 12 months from date of despatch |
| Acknowledgements | This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and |

non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](https://www.bio-rad.com/en-us/trademarks) for details.

His-tag is a registered trademark of EMD Biosciences.

Remicade is a trademark of Janssen Biotech, Inc.

| | |
|--------------------------------------|--|
| Health And Safety Information | Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA009P 20479 |
| Licensed Use | For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad. |
| Regulatory | For research purposes only |
| Technical Advice | Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual |

Related Products

Recommended Useful Reagents

[HUMAN ANTI INFLIXIMAB \(DRUG/TARGET COMPLEX\) \(TZA009\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA212\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA213\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA233\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA234\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA214\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA215\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA216\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M428785:240301'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)