

Datasheet: TZA0105

Description:	HUMAN ANTI IXEKIZUMAB
Specificity:	IXEKIZUMAB
Other names:	Taltz
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD60570ad
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL phage display library and expressed in a proprietary E. coli strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

Ixezumab

Specificity

Human Anti-Ixezumab Antibody, clone AbD60570ad is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug ixekizumab, and its biosimilar products. It does not recognize the drug target, human Interleukin 17 (IL-17A), nor ixekizumab in complex with IL-17A. This antibody can be used to measure free ixekizumab levels in serum from patients.

A pair of anti-ixekizumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD60570ad is recommended as a capture antibody, paired with Human Anti-Ixezumab Antibody, clone AbD60557 ([TZA0106](#)) as the detection antibody.

TZA0105 is a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the [SpyCatchers](#) available in our catalog.

Ixezumab (trade name Taltz) is a humanized monoclonal antibody of the IgG4/kappa isotype that binds to interleukin 17 (IL-17A), a pro-inflammatory cytokine, and blocks its action. Blocking its effect is beneficial for use in inflammatory conditions. Ixezumab is indicated for the treatment of moderate-to-severe plaque psoriasis and as a second-line therapy for active psoriatic arthritis.

[View a summary of all anti-ixekizumab antibodies](#)

Affinity

The monovalent intrinsic affinity of Human anti-Ixezumab Antibody, clone AbD60570ad was measured as $K_D = 0.2$ nM by real time, label-free molecular interaction analysis on immobilized ixekizumab.

ELISA

Clone AbD60570ad can be used in indirect ELISA to detect ixekizumab. It can also be used as a capture antibody to develop a pharmacokinetic (PK) bridging ELISA for risankizumab, together with [TZA0106](#) or [TZA0106P](#) as the detection reagent.

Protocol: [PK bridging ELISA](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Taltz is a registered trademark of Eli Lilly.

Health And Safety Information	Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA0105
Licensed Use	For <i>in vitro</i> . research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)
[HUMAN ANTI IXEKIZUMAB \(HCA432\)](#)
[HUMAN ANTI IXEKIZUMAB \(HCA433\)](#)
[HUMAN ANTI IXEKIZUMAB \(HCA434\)](#)
[HUMAN ANTI IXEKIZUMAB \(HCA435\)](#)
[HUMAN ANTI IXEKIZUMAB \(TZA0106\)](#)
[HUMAN ANTI IXEKIZUMAB:HRP \(TZA0106P\)](#)
[HUMAN ANTI IXEKIZUMAB \(TZA0107\)](#)
[HUMAN ANTI IXEKIZUMAB:HRP \(TZA0107P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M445594:250904'

Printed on 24 Sep 2025