

Datasheet: TZA0106P

**BATCH NUMBER 171364**

<b>Description:</b>	HUMAN ANTI IXEKIZUMAB:HRP
<b>Specificity:</b>	IXEKIZUMAB
<b>Other names:</b>	Taltz
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD60557
<b>Isotype:</b>	Fab antibody
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Product Form

A bivalent human recombinant Fab (kappa light chain) generated by coupling of the monovalent SpyTagged Fab [TZA0106](#) with a HRP conjugated BiSpyCatcher version 2. The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

### Preparation

Metal chelate affinity chromatography

### Source

E.coli

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.0095% MIT

### Approx. Protein

Antibody concentration 0.5 mg/ml

## Concentrations

---

**Immunogen** Ixekizumab

---

**Specificity** **Human Anti-Ixekizumab Antibody, clone AbD60557** is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug ixekizumab, and its biosimilar products. It does not recognize the drug target, human Interleukin 17 (IL-17A), nor ixekizumab in complex with IL-17A. This antibody can be used to measure free ixekizumab levels in serum from patients.

Clone AbD60557 is available as a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the [SpyCatchers](#) available in our catalog under product code [TZA0106](#). It is also available as a bivalent Fab antibody already conjugated to HRP via coupling to the HRP conjugated BiSpyCatcher2 ([TZC002P](#)) for use in direct ELISA under product code [TZA0106P](#).

A pair of anti-ixekizumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD60557 is recommended as a detection antibody, paired with Human Anti-Ixekizumab Antibody, clone AbD60570ad ([TZA0105](#)) as the capture reagent.

Ixekizumab (trade name Taltz) is a humanized monoclonal antibody of the IgG4/kappa isotype that binds to interleukin 17 (IL-17A), a pro-inflammatory cytokine, and blocks its action. Blocking its effect is beneficial for use in inflammatory conditions. Ixekizumab is indicated for the treatment of moderate-to-severe plaque psoriasis and as a second-line therapy for active psoriatic arthritis.

---

**Affinity** The monovalent intrinsic affinity of AbD60557pap was measured as  $K_D = 0.1$  nM by real time, label free molecular interaction analysis on immobilized ixekizumab.

---

**ELISA** Clone AbD60557pap can be used in direct ELISA to detect ixekizumab. It can also be used as a detection antibody to develop a pharmacokinetic (PK) bridging ELISA for ixekizumab, together with [TZA0105](#) as the capture reagent.

---

**Storage** This product is shipped frozen.  
When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.  
His-tag is a registered trademark of EMD Biosciences.  
Cynamza is a trademark of Eli Lilly.

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20479 available at: <a href="https://www.bio-rad-antibodies.com/SDS/TZA0106P">https://www.bio-rad-antibodies.com/SDS/TZA0106P</a>
<b>Licensed Use</b>	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only
<b>Technical Advice</b>	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at <a href="#">SpyTag and SpyCatcher Products</a> , and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a>

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)  
[HUMAN ANTI IXEKIZUMAB \(HCA432\)](#)  
[HUMAN ANTI IXEKIZUMAB \(HCA433\)](#)  
[HUMAN ANTI IXEKIZUMAB \(HCA434\)](#)  
[HUMAN ANTI IXEKIZUMAB \(HCA435\)](#)  
[HUMAN ANTI IXEKIZUMAB \(TZA0105\)](#)  
[HUMAN ANTI IXEKIZUMAB \(TZA0106\)](#)  
[HUMAN ANTI IXEKIZUMAB:HRP \(TZA0106P\)](#)  
[HUMAN ANTI IXEKIZUMAB \(TZA0107\)](#)  
[HUMAN ANTI IXEKIZUMAB:HRP \(TZA0107P\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M444932:250814'

Printed on 24 Sep 2025