

Datasheet: TZA0123P

Description:	HUMAN ANTI DULAGLUTIDE
Specificity:	Dulaglutide
Other names:	Trulicity
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD64818
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A bivalent human recombinant Fab (lambda light chain) generated by coupling of the monovalent SpyTagged Fab [TZA0123](#) with a HRP conjugated BiSpyCatcher version 2. The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen	Dulaglutide
Specificity	<p>Human Anti-Dulaglutide Antibody, clone AbD64818 is specific for human GLP-1 and its peptidase-IV-protected analog present in the human IgG4-Fc fusion protein dulaglutide.</p> <p>Clone AbD64818 is available as a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the SpyCatchers available in our catalog under product code TZA0123. It is also available as a bivalent Fab antibody already conjugated to HRP via coupling to the HRP conjugated BiSpyCatcher2 (TZC002P) for use in direct ELISA under product code TZA0123P.</p> <p>A pair of anti-dulaglutide antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD64818ad is recommended as a capture antibody, paired with Human Anti-Dulaglutide Antibody, clone AbD64818pap (TZA0123P) as the detection reagent.</p> <p>Dulaglutide (trade name Trulicity) is a fusion protein, combining a modified GLP-1 analogue with a human IgG4-Fc heavy chain. It binds to GLP-1 and activates the GLP-1 receptor found in the pancreatic beta cells increasing intracellular cyclic AMP in beta cells, leading to insulin release and subsequent reduction of blood glucose concentrations. Dulaglutide is indicated for the management of type 2 diabetes and reduce major adverse cardiovascular events.</p> <p>View a summary of all anti-dulaglutide antibodies</p>
Affinity	The monovalent intrinsic affinity of AbD64818pap was measured as $K_D = 2.0$ nM by real time, label free molecular interaction analysis on immobilized dulaglutide.
ELISA	<p>Clone AbD64818pap can be used as a detection antibody to develop a pharmacokinetic (PK) bridging ELISA for dulaglutide, together with TZA0123 as the capture reagent.</p> <p>Protocol: PK bridging ELISA</p>
Storage	<p>This product is shipped frozen.</p> <p>When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p> <p>Trulicity is a trademark of Eli Lilly.</p>
Health And Safety	Material Safety Datasheet documentation #20479 available at:

Information <https://www.bio-rad-antibodies.com/SDS/TZA0123P>

Licensed Use For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at [SpyTag and SpyCatcher Products](#), and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI DULAGLUTIDE \(TZA0123\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M450609:260401'

Printed on 19 May 2026