

## Datasheet: TZA0125P

**BATCH NUMBER 172591**

<b>Description:</b>	HUMAN ANTI ISATUXIMAB:HRP
<b>Specificity:</b>	ISATUXIMAB
<b>Other names:</b>	Sarclisa
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD64592
<b>Isotype:</b>	Fab antibody
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Product Form

A bivalent human recombinant Fab (lambda light chain) generated by coupling of the monovalent SpyTagged Fab [TZA0125](#) with a HRP conjugated BiSpyCatcher version 2 ([TZC002P](#)). The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

### Preparation

Metal chelate affinity chromatography

### Source

E.coli

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.0095% MIT

### Approx. Protein

Antibody concentration 0.5 mg/ml

## Concentrations

---

**Immunogen** Isatuximab

---

**Specificity** **Human Anti-Isatuximab Antibody, clone AbD64592** is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug isatuximab. It does not recognize the drug target, CD38, nor isatuximab in complex with CD38. This antibody can be used to measure free isatuximab levels in serum from patients.

Clone AbD64592 is available as a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the [SpyCatchers](#) available in our catalog under product code [TZA0125](#). It is also available as a bivalent Fab antibody already conjugated to HRP via coupling to the HRP conjugated BiSpyCatcher2 ([TZC002P](#)) for use in direct ELISA under product code [TZA0125P](#).

A pair of anti-isatuximab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD64592 is recommended as a detection antibody, paired with Human Anti-Isatuximab Antibody, clone AbD64599ad ([TZA0124](#)) as the capture reagent.

Isatuximab (trade name Sarclisa) is a chimeric human/mouse IgG1/kappa, produced using recombinant DNA technology in Chinese hamster ovary (CHO) cells. Isatuximab is an IgG1-derived monoclonal antibody targeted against CD38. Its binding towards CD38 positive myeloma cells results in several downstream effects, including direct apoptosis of the affected cell and activation of immune mechanisms including antibody-dependent cell-mediated cytotoxicity (ADCC), antibody-dependent cellular phagocytosis (ADCP), and complement dependent cyto-toxicity (CDC), all of which result in potent anti-tumour activity. Isatuximab is indicated for the treatment of multiple myeloma in adults.

[View a summary of all anti-isatuximab antibodies](#)

---

**Affinity** The monovalent intrinsic affinity of Human Anti-Isatuximab Antibody, clone AbD64592pap was measured as  $K_D = 1.0$  nM by real time, label-free molecular interaction analysis on immobilized isatuximab.

---

**ELISA** Clone AbD64592pap can be used in direct ELISA to detect isatuximab. It can also be used as a detection antibody to develop a pharmacokinetic (PK) bridging ELISA for isatuximab, together with [TZA0124](#) as the capture reagent.

Protocol: [PK bridging ELISA](#)

---

**Storage** This product is shipped frozen.  
When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.

---

<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See <a href="http://bio-rad.com/en-us/trademarks">bio-rad.com/en-us/trademarks</a> for details.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p> <p>Sarclisa is a registered trademark of ImmunoGen, Sanofi, and Sanofi-Aventis U.S.LLC.</p>
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20479 available at: <a href="https://www.bio-rad-antibodies.com/SDS/TZA0125P">https://www.bio-rad-antibodies.com/SDS/TZA0125P</a>
<b>Licensed Use</b>	For <i>in vitro</i> . research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only
<b>Technical Advice</b>	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at <a href="#">SpyTag and SpyCatcher Products</a> , and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a> .

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI ISATUXIMAB \(TZA0124\)](#)

[HUMAN ANTI ISATUXIMAB \(TZA0125\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M450383:260318'

**Printed on 22 May 2026**