

Datasheet: TZA0126

Description:	HUMAN ANTI TREMELIMUMAB
Specificity:	TREMELIMUMAB
Other names:	Imjudo
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD64576ad
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in a proprietary *E. coli* strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E. coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

Tremelimumab

Specificity

Human Anti-Tremelimumab Antibody, clone AbD64576ad is a paratope specific, inhibitory anti-idiotypic antibody (Type 1) that specifically recognizes the humanized monoclonal antibody drug tremelimumab. It does not recognize the drug target CTLA-4, nor tremelimumab in complex with CTLA-4. This antibody can be used to measure free tremelimumab levels in serum from patients.

A pair of anti-tremelimumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD64576ad is recommended as a capture antibody, paired with Human Anti-Tremelimumab Antibody, clone AbD64580pap ([TZA0127P](#)) as the detection antibody.

TZA0126 is a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the [SpyCatchers](#) available in our catalog.

Tremelimumab (trade name Imjudo) is a fully humanized monoclonal IgG2/kappa antibody with allotype IGHG2*01. It binds to cytotoxic T-lymphocyte-associated protein 4 (CTLA-4) to enhance T cell-mediated killing of tumours and reduces tumour growth. Tremelimumab in combination with durvalumab and platinum-based chemo-therapy is indicated for the first-line treatment of adults with metastatic non-small cell lung cancer (NSCLC) with no sensitising EGFR mutations or ALK positive mutations.

[View a summary of all anti-tremelimumab antibodies](#)

Affinity

The monovalent intrinsic affinity of AbD64576ad was measured as $K_D = < 0.02$ nM by real time, label free molecular interaction analysis on immobilized tremelimumab.

ELISA

Clone AbD64576ad can be used in indirect ELISA to detect tremelimumab. It can also be used as a capture antibody to develop a pharmacokinetic (PK) bridging ELISA for tremelimumab together with [TZA0127](#) or [TZA0127P](#) as the detection reagent.

Protocol: [PK bridging ELISA](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

His-tag is a registered trademark of EMD Biosciences.

Imjudo is a trademark of Pfizer and now investigated by MedImmune (AstraZeneca).

Health And Safety Information	Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA0126
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI TREMELIMUMAB \(TZA0127\)](#)

[HUMAN ANTI TREMELIMUMAB:HRP \(TZA0127P\)](#)

[HUMAN ANTI TREMELIMUMAB \(DRUG/TARGET COMPLEX\) \(TZA0128\)](#)

[HUMAN ANTI TREMELIMUMAB \(DRUG/TARGET COMPLEX\):HRP \(TZA0128P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M450595:260330'

Printed on 30 Mar 2026