

Datasheet: TZA0137P

Description:	HUMAN ANTI YTE MUTATIONS
Specificity:	YTE MUTATIONS
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD61602
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A bivalent human recombinant Fab (lambda light chain) generated by coupling of the monovalent SpyTagged Fab [TZA0137](#) with a HRP conjugated BiSpyCatcher version 2. The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

YTE mutations (AbD18705iu, hIgG1 carrying the YTE mutations)

Specificity	<p>Human Anti-YTE mutation Antibody, clone AbD61602 is an antibody that specifically recognizes the YTE mutation (M252Y/S254T/T256E) in human IgG1 antibodies. This antibody can be used to measure antibody drug concentrations in serum from patients treated with antibody drugs carrying this YTE mutation.</p> <p>For PK assays, clone AbD61602ad is recommended as a capture antibody, paired with Human anti-YTE mutation antibody, clone AbD61602pap (TZA0137P) as the detection antibody. It can also be paired with anti-idiotypic detection antibodies of the drug to be measured.</p> <p>Clone AbD61602 is available as a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the SpyCatchers available in our catalog under product code TZA0137. It is also available as a bivalent Fab antibody already conjugated to HRP via coupling to the HRP conjugated BiSpyCatcher2 (TZC002P) for use in direct ELISA under product code TZA0137P.</p> <p>The YTE-mutations (M252Y/S254T/T256E) were initially introduced at the CH2-CH3 interface of the humanized IgG1 drug antibody MEDI-524 (anti-respiratory syncytial virus, RSV). It was shown that these mutations lead to a prolonged half time in serum due to a higher affinity to the neonatal Fc receptor (FcRn) and is therefore an interesting format for many antibody drugs.</p> <p>View a summary of all anti-yte mutation antibodies</p>
Affinity	<p>The monovalent intrinsic affinity of AbD61602pap was measured as $K_D = 7$ nM by real time, label free molecular interaction analysis on immobilized AbD18705iu.</p>
ELISA	<p>Clone AbD61602pap can be used in indirect ELISA for AbD18705iu and other antibodies carrying the YTE mutations. It can also be used as a detection antibody to develop a pharmacokinetic (PK) bridging ELISA together with TZA0137 as the capture reagent. Note: This antibody is recommended to use when working with antibodies carrying other Fc mutations in addition to the YTE mutations.</p> <p>Protocol: PK bridging ELISA</p>
Storage	<p>This product is shipped frozen.</p> <p>When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.</p>
Guarantee	<p>12 months from date of despatch</p>
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p>
Health And Safety	<p>Material Safety Datasheet documentation #20479 available at:</p>

Information	https://www.bio-rad-antibodies.com/SDS/TZA0137P
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI YTE MUTATIONS \(TZA0135\)](#)

[HUMAN ANTI YTE MUTATIONS \(TZA0136\)](#)

[HUMAN ANTI YTE MUTATIONS \(TZA0136P\)](#)

[HUMAN ANTI YTE MUTATIONS \(TZA0137\)](#)

[RECOMBINANT HUMAN IgG1 YTE KAPPA ALLOTYPIC G1m3 \(HCA413\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M450664:260420'

Printed on 14 May 2026