

Datasheet: TZA0145P

BATCH NUMBER 172701

Description:	HUMAN ANTI ETANERCEPT:HRP
Specificity:	ETANERCEPT
Other names:	Enbrel
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD26183pap
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A bivalent human recombinant Fab (kappa light chain) generated by coupling of the monovalent SpyTagged Fab with a HRP conjugated BiSpyCatcher version 2 ([TZC002P](#)). The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein

Antibody concentration 0.5 mg/ml

Concentrations

Immunogen Etanercept

Specificity **Human Anti-Etanercept Antibody, clone AbD26183pap** is a non-inhibitory anti-idiotypic antibody (type 2) that specifically recognizes free humanized monoclonal antibody drug etanercept. It does not recognize free TNF- α and TNF- β but Etanercept in complex with TNF- α . This antibody can be used to measure free and bound Etanercept in serum of patients.

TZA0145P, clone AbD26183pap is a bivalent Fab antibody already conjugated to HRP via coupling to the HRP conjugated BiSpyCatcher2 ([TZC002P](#)) for use in direct ELISA.

A pair of anti-etanercept antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. Clone AbD26183pap is recommended as the detection antibody, paired with Human Anti-Etanercept Antibody in monovalent Fab format, clone AbD25939 ([HCA276](#)) as the capture antibody.

Etanercept (Enbrel[®]) is a dimeric fusion protein consisting of the extracellular ligand-binding portion of the human 75 kDa (p75) TNFR2 linked to the Fc portion of human IgG1. The Fc component of etanercept contains the CH2 domain, the CH3 domain and hinge region, but not the CH1 domain of IgG1. It is used to treat autoimmune diseases by interfering with TNF, acting as a TNF inhibitor. It binds to both TNF α and TNF β (Lymphotoxin A).

[View a summary of all anti-etanercept antibodies.](#)

Affinity The intrinsic affinity of the monovalent form of Human anti Etanercept antibody, clone AbD26183pap is $K_D = 14.8$ nM as measured by real time, label free molecular interaction analysis on immobilized Etanercept.

ELISA Clone AbD26183pap can be used in direct or indirect ELISA to detect etanercept. It can also be used as a detection antibody to develop a pharmacokinetic (PK) bridging ELISA for etanercept together with [HCA276](#) as the capture reagent.

Protocol: [PK bridging ELISA](#)

Storage This product is shipped frozen.
When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week

Guarantee 12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Enbrel is a registered trademark of Immunex Corporation.

Health And Safety Information	Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA0145P
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)
[HUMAN ANTI ETANERCEPT \(HCA277\)](#)
[HUMAN ANTI ETANERCEPT \(HCA276\)](#)
[HUMAN ANTI ETANERCEPT:HRP \(HCA277P\)](#)
[HUMAN ANTI ETANERCEPT \(HCA279\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M450066:260303'

Printed on 19 Mar 2026