

## Datasheet: TZA0147

<b>Description:</b>	HUMAN ANTI NIRSEVIMAB
<b>Specificity:</b>	NIRSEVIMAB
<b>Other names:</b>	Beyfortus
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD64623ad
<b>Isotype:</b>	Fab antibody
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in a proprietary *E. coli* strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

### Preparation

Metal chelate affinity chromatography

### Source

*E. coli*

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.0095% MIT

### Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

<b>Immunogen</b>	Nirsevimab
<b>Specificity</b>	<p><b>Human Anti-Nirsevimab Antibody, clone AbD64623ad</b> is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug nirsevimab, and its biosimilar products. It does not recognize the free Fusion glycoprotein F0 of RSV, nor nirsevimab in complex with Fusion glycoprotein F0 of RSV. This antibody can be used to measure free nirsevimab levels in serum from patients.</p> <p>A pair of anti-nirsevimab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD64623ad is recommended as a capture antibody, paired with Human Anti-Nirsevimab Antibody, clone AbD64624pap (<a href="#">TZA0148P</a>) as the detection antibody.</p> <p>TZA0147 is a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the <a href="#">SpyCatchers</a> available in our catalog.</p> <p>Nirsevimab (trade name Beyfortus) is a humanized monoclonal antibody of the IgG1/kappa isotype that targets the prefusion conformation of the RSV F protein, a glycoprotein involved in the membrane fusion step of the viral entry process and neutralises several RSV A and B strains. Nirsevimab is indicated to treat respiratory syncytial virus (RSV) lower respiratory tract disease in neonates and infants.</p>
<b>Affinity</b>	The monovalent intrinsic affinity of AbD64623ad was measured as $K_D = 0.1$ nM by real time, label free molecular interaction analysis on immobilized nirsevimab.
<b>ELISA</b>	Clone AbD64623ad can be used in indirect ELISA to detect nirsevimab. It can also be used as a capture antibody to develop a pharmacokinetic (PK) bridging ELISA for nirsevimab together with <a href="#">TZA0148</a> or <a href="#">TZA0148P</a> as the detection reagent.
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See <a href="http://bio-rad.com/en-us/trademarks">bio-rad.com/en-us/trademarks</a> for details.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p> <p>Beyfortus is a trademark of Sanofi and AstraZeneca.</p>
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20479 available at: <a href="https://www.bio-rad-antibodies.com/SDS/TZA0147">https://www.bio-rad-antibodies.com/SDS/TZA0147</a>

<b>Licensed Use</b>	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only
<b>Technical Advice</b>	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at <a href="#">SpyTag and SpyCatcher Products</a> , and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a>

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI NIRSEVIMAB \(TZA0148\)](#)

[HUMAN ANTI NIRSEVIMAB:HRP \(TZA0148P\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M449475:260126'

**Printed on 29 Jan 2026**