

Datasheet: TZA025

Description:	HUMAN ANTI PERTUZUMAB
Specificity:	PERTUZUMAB
Other names:	Perjeta
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD35569
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (kappa light chain) selected from the HuCAL® phage display library, expressed in a proprietary *E. coli* strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E. coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

Pertuzumab

Specificity

Human Anti-Pertuzumab Antibody, clone AbD35569 is a non-paratope specific, non-inhibitory anti-idiotypic antibody that specifically recognizes the humanized monoclonal antibody drug pertuzumab. It recognizes free pertuzumab and pertuzumab in complex with its target, the human epidermal growth factor receptor 2 (HER2) but not free HER2. It can be used to measure total levels of pertuzumab and biosimilars in patient samples.

The monovalent Fab format of this antibody, clone AbD35569ad ([TZA025](#)) is tagged with a SpyTag2 at the C-terminus of the Fab heavy chain, enabling the user to couple it to a [SpyCatcher reagent](#) for conversion to alternative formats in less than one hour.

This antibody is also available as a bivalent Fab antibody pre-coupled to the HRP-conjugated BiSpyCatcher2 ([TZC002P](#)) under product code [TZA025P](#) for direct detection of pertuzumab in bioanalysis assays.

A pair of anti-pertuzumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure total drug. Clone AbD35569 can be used as a detection antibody to develop a pharmacokinetic (PK) bridging assay together with [HCA341](#) as the capture antibody.

Pertuzumab (Perjeta) is an IgG1/kappa antibody humanized from mouse. It has been approved for the treatment of patients with metastatic HER2-positive breast cancer in combination with trastuzumab. Pertuzumab binds to the extracellular domain II of HER2 and inhibits HER2–HER3 dimerization, reducing signaling through intracellular pathways ([Capelan et al. 2013](#)).

[View a summary of all anti-pertuzumab antibodies](#)

Affinity

The monovalent intrinsic affinity of AbD35569ad was measured as $K_D = 31$ nM by real time, label free molecular interaction analysis on immobilized pertuzumab.

ELISA

Clone AbD35569ad can be used in indirect ELISA to detect pertuzumab. It can also be used as a detection antibody, together with [HCA341](#) as the capture antibody to develop a PK bridging ELISA to detect pertuzumab total drug levels in the serum of patients.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

His-tag is a registered trademark of EMD Biosciences.
Perjeta is a trademark of Genentech, Inc.

Health And Safety Information	Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA025 20479
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[SpyTag3 PEPTIDE \(BLP086\)](#)
[HISPEC ASSAY DILUENT \(BUF049A\)](#)
[BiSpyCatcher2:HRP \(TZC002P\)](#)
[HUMAN ANTI PERTUZUMAB \(HCA341\)](#)
[HUMAN ANTI PERTUZUMAB \(HCA342\)](#)
[HUMAN ANTI PERTUZUMAB \(HCA344\)](#)
[HUMAN ANTI PERTUZUMAB \(HCA345\)](#)
[HUMAN ANTI PERTUZUMAB:HRP \(TZA025P\)](#)
[HUMAN ANTI PERTUZUMAB \(DRUG/TARGET COMPLEX\) \(TZA035\)](#)
[HUMAN ANTI PERTUZUMAB \(DRUG/TARGET COMPLEX\):HRP \(TZA035P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M432866:240912'

Printed on 12 Sep 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)