

Datasheet: TZA025P

BATCH NUMBER 168431

Description:	HUMAN ANTI PERTUZUMAB:HRP
Specificity:	PERTUZUMAB
Other names:	Perjeta
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD35569
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A bivalent human recombinant Fab (kappa light chain) generated by coupling of the monovalent SpyTagged Fab [TZA025](#) with a HRP conjugated BiSpyCatcher version 2 ([TZC002P](#)). The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein

Antibody concentration 0.5 mg/ml

Concentrations

Immunogen Pertuzumab

Specificity **Human Anti-Pertuzumab Antibody, clone AbD35569** is a non-paratope specific, non-inhibitory anti-idiotypic antibody that specifically recognizes the humanized monoclonal antibody drug pertuzumab. It recognizes free pertuzumab and pertuzumab in complex with its target, the human epidermal growth factor receptor 2 (HER2) but not free HER2. It can be used to measure total levels of pertuzumab and biosimilars in patient samples.

The monovalent Fab format of this antibody, clone AbD35569ad ([TZA025](#)) is tagged with a SpyTag2 at the C-terminus of the Fab heavy chain, enabling the user to couple it to a [SpyCatcher reagent](#) for conversion to alternative formats in less than one hour.

This antibody is also available as a bivalent Fab antibody pre-coupled to the HRP-conjugated BiSpyCatcher2 ([TZC002P](#)) under product code [TZA025P](#) for direct detection of pertuzumab in bioanalysis assays.

A pair of anti-pertuzumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure total drug. Clone AbD35569 can be used as a detection antibody to develop a pharmacokinetic (PK) bridging assay together with [HCA341](#) as the capture antibody.

Pertuzumab (Perjeta) is an IgG1/kappa antibody humanized from mouse. It has been approved for the treatment of patients with metastatic HER2-positive breast cancer in combination with trastuzumab. Pertuzumab binds to the extracellular domain II of HER2 and inhibits HER2–HER3 dimerization, reducing signaling through intracellular pathways ([Capelan *et al.* 2013](#)).

[View a summary of all anti-pertuzumab antibodies](#)

Affinity The monovalent intrinsic affinity of AbD35569pap was measured as $K_D = 31$ nM by real time, label free molecular interaction analysis on immobilized pertuzumab.

ELISA Clone AbD35569pap can be used in direct ELISA to detect pertuzumab. It can also be used as a detection antibody, together with [HCA341](#) as the capture antibody to develop a PK bridging ELISA to detect pertuzumab total drug levels in the serum of patients.

Storage This product is shipped frozen.
When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.

Guarantee 12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

His-tag is a registered trademark of EMD Biosciences.

Perjeta is a trademark of Genentech, Inc.

Health And Safety Information	Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA025P
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA341\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA342\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA344\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA345\)](#)

[HUMAN ANTI PERTUZUMAB \(TZA025\)](#)

[HUMAN ANTI PERTUZUMAB \(DRUG/TARGET COMPLEX\) \(TZA035\)](#)

[HUMAN ANTI PERTUZUMAB \(DRUG/TARGET COMPLEX\):HRP \(TZA035P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M431717:240813'

Printed on 20 Jun 2025