

Datasheet: TZA027P

BATCH NUMBER 164327

Description:	HUMAN ANTI HUMAN IgD:HRP
Specificity:	IgD
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD39958
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

A bivalent human recombinant Fab (lambda light chain) generated by coupling of the monovalent SpyTagged Fab [TZA027](#) with a HRP conjugated BiSpyCatcher version 2. The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Native human IgD
Specificity	<p>Human anti human IgD, clone AbD39958 recognizes human IgD. No cross-reactivity is seen with other human immunoglobulin subclasses.</p> <p>The monovalent Fab format of this antibody (TZA027) is tagged with a SpyTag at the C-terminus of the Fab heavy chain, enabling the user to couple it to a SpyCatcher reagent for conversion to alternative formats in less than an hour.</p> <p>Human anti human IgD, clone AbD39958 is also available as a bivalent Fab antibody already coupled to HRP (TZA027P).</p>
Storage	<p>This product is shipped frozen.</p> <p>When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details</p> <p>His-tag is a registered trademark of EMD Biosciences</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA027P</p> <p>20479</p>
Regulatory	For research purposes only
Technical Advice	<p>Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products, and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual</p>

Related Products

Recommended Useful Reagents

[SpyTag3 PEPTIDE \(BLP086\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M418499:230424'

Printed on 01 Mar 2024

