

Datasheet: TZA029

BATCH NUMBER 164332

Description:	HUMAN ANTI HUMAN IgM
Specificity:	IgM
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD39983
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in *E. coli* strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Recombinant human IgM
Specificity	<p>Human anti human IgM, clone AbD39983 recognizes human IgM. No cross-reactivity is seen with other human immunoglobulin subclasses.</p> <p>The monovalent Fab format of this antibody (TZA029) is tagged with a SpyTag at the C-terminus of the Fab heavy chain, enabling the user to couple it to a SpyCatcher reagent for conversion to alternative formats in less than an hour.</p> <p>Human anti human IgM, clone AbD39983 is also available as a bivalent Fab antibody already coupled to HRP (TZA029P).</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details</p> <p>His-tag is a registered trademark of EMD Biosciences</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA029</p> <p>20479</p>
Regulatory	For research purposes only
Technical Advice	<p>Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products, and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual</p>

Related Products

Recommended Useful Reagents

[SpyTag3 PEPTIDE \(BLP086\)](#)

North & South America
 Tel: +1 800 265 7376
 Fax: +1 919 878 3751
 Email: antibody_sales_us@bio-rad.com

Worldwide
 Tel: +44 (0)1865 852 700
 Fax: +44 (0)1865 852 739
 Email: antibody_sales_uk@bio-rad.com

Europe
 Tel: +49 (0) 89 8090 95 21
 Fax: +49 (0) 89 8090 95 50
 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
 'M418501:230424'

