

Datasheet: TZA036P

Description:	HUMAN ANTI RAT IgG1:HRP
Specificity:	IgG1
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD41031
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	A bivalent human recombinant Fab (lambda light chain) generated by coupling of the monovalent SpyTagged Fab TZA036 with a HRP conjugated BiSpyCatcher version 2. The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide (BLP086) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.
Preparation	Metal chelate affinity chromatography
Source	E.coli
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.0095% MIT
Approx. Protein	Antibody concentration 0.5 mg/ml

Concentrations

Immunogen	Recombinant human/rat IgG1 chimera
------------------	------------------------------------

External Database Links

UniProt:

[P20759](#)

[Related reagents](#)

Entrez Gene:

[299354](#)

Ighg

[Related reagents](#)

Specificity

Human anti rat IgG1 antibody, clone AbD41031 recognizes rat IgG1. No cross-reactivity is seen with other rat immunoglobulin subclasses.

The monovalent Fab format of this antibody (TZA036) is tagged with a SpyTag at the C-terminus of the Fab heavy chain, enabling the user to couple it to a [SpyCatcher reagent](#) for conversion to alternative formats in less than an hour.

Human anti rat IgG1, clone AbD41031 is also available as a bivalent Fab antibody already coupled to HRP (TZA036P).

Storage

This product is shipped frozen.

When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details
His-tag is a registered trademark of EMD Biosciences

Health And Safety Information

Material Safety Datasheet documentation #20479 available at:
<https://www.bio-rad-antibodies.com/SDS/TZA036P>

Regulatory

For research purposes only

Technical Advice

Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at [SpyTag and SpyCatcher Products](#), and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M450163:260306'

Printed on 06 Mar 2026

