

## Datasheet: TZA044

<b>Description:</b>	HUMAN ANTI OCRELIZUMAB
<b>Specificity:</b>	OCRELIZUMAB
<b>Other names:</b>	Ocrevus
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD50641ad
<b>Isotype:</b>	Fab antibody
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Product Form

A monovalent human recombinant Fab (lambda light chain) selected from Bio-Rad's human antibody phage display libraries. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

### Preparation

Metal chelate affinity chromatography

### Source

E.coli

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.0095% MIT

### Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

**Specificity**

**Human Anti-Ocrelizumab Antibody, clone AbD50641ad** is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug ocrelizumab and its biosimilar products. Consequently, this antibody can be used to measure free ocrelizumab levels in serum from patients.

This antibody is available in two formats - a monovalent Fab unconjugated format, clone AbD50641ad (TZA044), which has a SpyTag incorporated. This format can be used as a capture antibody in a pharmacokinetic (PK) bridging ELISA together with AbD50651 as a detection antibody. Clone AbD50641 is also available in full length IgG1 format under product code [HCA400](#) which can be used in an ADA bridging ELISA as a surrogate positive control.

Human Anti-Ocrelizumab Antibody, clone AbD50641ad has been shown to cross-react with the antibody drug rituximab (trade name Rituxan). Ocrelizumab is the humanized form of rituximab. There are changes in at least one of the CDRs and it contains Fc modification.

Ocrelizumab (trade name Ocrevus) is a monoclonal IgG1/kappa antibody (VH3/ $\kappa$ 1). In vitro, it binds to CD20 cell surface molecules on target cells thereby inducing immune effector mechanisms such as antibody-dependent cellular phagocytosis (ADCP), antibody-dependent cellular cytotoxicity (ADCC), complement-dependent cytotoxicity (CDC), and apoptosis. In vivo, ocrelizumab selectively and effectively depletes CD20<sup>+</sup> B cells. Hence, it is an immunosuppressive drug. Ocrelizumab is used to treat primary progressive multiple sclerosis and relapsing forms of multiple sclerosis in adults (including clinically isolated syndrome, relapsing-remitting disease, and active secondary progressive disease).

[View a summary of all-ocrelizumab antibodies](#)

**Affinity**

The monovalent intrinsic affinity of Human anti-Ocrelizumab Antibody, clone AbD50641ad was measured as  $K_D < 0.001$  nM by real time, label-free molecular interaction analysis on immobilized ocrelizumab.

**ELISA**

This product may be used in a direct ELISA.  
Clone AbD50641ad can be used as the capture antibody for ocrelizumab in a PK bridging ELISA together with clone AbD50651pap ([TZA045P](#)) as the detection antibody.

Protocol: [PK bridging ELISA](#)

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See <a href="http://bio-rad.com/en-us/trademarks">bio-rad.com/en-us/trademarks</a> for details.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p> <p>Ocrevus is a registered trademark of Genentech/ Roche and Biogen.</p> <p>Rituxan is a registered trademark of Genentech and Biogen.</p>
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #20479 available at: <a href="https://www.bio-rad-antibodies.com/SDS/TZA044">https://www.bio-rad-antibodies.com/SDS/TZA044</a></p> <p>20479</p>
<b>Licensed Use</b>	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only
<b>Technical Advice</b>	<p>Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at <a href="#">SpyTag and SpyCatcher Products</a>, and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a></p>

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)  
[SpyTag3 PEPTIDE \(BLP086\)](#)  
[HUMAN ANTI OCRELIZUMAB \(HCA400\)](#)  
[HUMAN ANTI OCRELIZUMAB \(HCA401\)](#)  
[HUMAN ANTI OCRELIZUMAB \(TZA045\)](#)  
[HUMAN ANTI OCRELIZUMAB:HRP \(TZA045P\)](#)  
[BiSpyCatcher2:HRP \(TZC002P\)](#)  
[HUMAN ANTI OCRELIZUMAB \(HCA399\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M442142:250528'

Printed on 17 Jun 2025