

Datasheet: TZA071

Description:	HUMAN ANTI GUSELKUMAB
Specificity:	GUSELKUMAB
Other names:	Tremfya
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD59937
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in a proprietary *E. coli* strain. The antibody is tagged with a DYKDDDDK tag, a SpyTag version 2 (VPTIVMVDAYKRYK) and a His6-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Protocol: [PK bridging ELISA](#)

Preparation	Metal chelate affinity chromatography
Source	<i>E.coli</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.0095% MIT

Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Guselkumab
Specificity	<p>Human Anti-Guselkumab Antibody, clone AbD59937 is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug guselkumab, and its biosimilar products. It does not recognize the drug target, the human Interleukin-23 (IL-23), nor guselkumab in complex with IL-23. This antibody can be used to measure free guselkumab levels in serum from patients.</p> <p>A pair of anti-guselkumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug; clone AbD59937 is recommended as the detection reagent, paired with Human Anti-Guselkumab Antibody, clone AbD59307ad (TZA070) as the capture antibody.</p> <p>TZA071 is a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the SpyCatchers available in our catalog.</p> <p>Clone AbD59937 is also available conjugated to HRP via coupling to the BiSpyCatcher2:HRP (TZC002P) for direct detection of guselkumab.</p> <p>Guselkumab (trade name Tremfya) is a fully human monoclonal antibody of the IgG1/lambda isotype binding to the p19 subunit of IL-23 to block its interaction with the IL-23 receptor, subsequently inhibiting the release of pro-inflammatory cytokines and chemokines. Guselkumab is approved for the treatment of moderate to severe plaque psoriasis in adults.</p> <p>View a summary of all anti-guselkumab antibodies</p>
Affinity	The monovalent intrinsic affinity of AbD59937ad was measured as $K_D = 0.06$ nM by real time, label free molecular interaction analysis on immobilized guselkumab.
ELISA	Clone AbD59937ad can be used in indirect ELISA to detect guselkumab. It can also be used as a detection antibody to develop a pharmacokinetic (PK) bridging ELISA for guselkumab, together with TZA070 as the capture reagent.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See

bio-rad.com/en-us/trademarks for details.
His-tag is a registered trademark of EMD Biosciences.
Tremfya is a trademark of MorphoSys.

Health And Safety Information	Material Safety Datasheet documentation #20479 available at: https://www.bio-rad-antibodies.com/SDS/TZA071
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at SpyTag and SpyCatcher Products , and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI GUSELKUMAB \(HCA420\)](#)

[HUMAN ANTI GUSELKUMAB \(HCA421\)](#)

[HUMAN ANTI GUSELKUMAB \(TZA070\)](#)

[HUMAN ANTI GUSELKUMAB:HRP \(TZA071P\)](#)

[HUMAN ANTI GUSELKUMAB \(DRUG/TARGET COMPLEX\) \(TZA072\)](#)

[HUMAN ANTI GUSELKUMAB \(DRUG/TARGET COMPLEX\):HRP \(TZA072P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M444669:250717'

Printed on 24 Sep 2025