

Datasheet: TZA072P

BATCH NUMBER 169961

Description:	HUMAN ANTI GUSELKUMAB (DRUG/TARGET COMPLEX):HRP
Specificity:	GUSELKUMAB DRUG/TARGET COMPLEX
Other names:	Tremfya
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD59967
Isotype:	Fab antibody
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A bivalent human recombinant Fab (kappa light chain) generated by coupling of the monovalent SpyTagged Fab [TZA072](#) with a HRP conjugated BiSpyCatcher version 2 ([TZC002P](#)). The coupled antibody has two DYKDDDDK tags and three His6-tags. It contains SpyTag3 peptide ([BLP086](#)) in 5-fold molar excess to block any unreacted Catcher sites. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Source

E.coli

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.0095% MIT

Approx. Protein

Antibody concentration 0.5 mg/ml

Concentrations

Immunogen Guselkumab

Specificity **Human anti guselkumab (drug/target complex) Antibody, clone AbD59967** recognizes the antibody drug guselkumab when in complex with its target, the human Interleukin-23 (IL-23). It does not recognize free guselkumab nor unbound IL-23.

Clone AbD59967 can be used to develop a pharmacokinetic (PK) antigen capture assay to measure levels of guselkumab captured via immobilized IL-23.

TZA072 is a monovalent Fab antibody incorporating a SpyTag at the C-terminus end of the antibody heavy chain for conversion by the end user into alternative formats using any of the [SpyCatchers](#) available in our catalog.

Clone AbD59967 is also available conjugated to HRP via coupling to the BiSpyCatcher2:HRP ([TZC002P](#)) for direct detection of the IL-23/guselkumab complex.

Guselkumab (trade name Tremfya) is a fully human monoclonal antibody of the IgG1/lambda isotype binding to the p19 subunit of IL-23 to block its interaction with the IL-23 receptor, subsequently inhibiting the release of pro-inflammatory cytokines and chemokines. Guselkumab is approved for the treatment of moderate to severe plaque psoriasis in adults.

Affinity The monovalent intrinsic affinity of AbD59967pap was measured as $K_D = 33$ nM by real time, label free molecular interaction analysis on immobilized IL-23A/guselkumab complex.

ELISA Clone AbD59967pap can be used in direct ELISA to detect the IL-23A/ guselkumab complex in a pharmacokinetic (PK) antigen capture assay.

Storage This product is shipped frozen.
When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.

Guarantee 12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

His-tag is a registered trademark of EMD Biosciences.

Tremfya is a trademark of MorphoSys.

Health And Safety Information Material Safety Datasheet documentation #20479 available at: <https://www.bio-rad-antibodies.com/SDS/TZA072P>

Licensed Use For *in vitro* research purposes and for commercial applications for the provision of *in vitro* testing services to support preclinical and clinical drug development. Any re-sale in any

form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols for coupling a Fab antibody with a SpyTag to a SpyCatcher can be found at [SpyTag and SpyCatcher Products](#), and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI GUSELKUMAB \(TZA071\)](#)

[HUMAN ANTI GUSELKUMAB:HRP \(TZA071P\)](#)

[HUMAN ANTI GUSELKUMAB \(TZA070\)](#)

[HUMAN ANTI GUSELKUMAB \(DRUG/TARGET COMPLEX\) \(TZA072\)](#)

[HUMAN ANTI GUSELKUMAB \(HCA420\)](#)

[HUMAN ANTI GUSELKUMAB \(HCA421\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M434357:250217'

Printed on 24 Sep 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)