

Datasheet: TZA077GA

BATCH NUMBER 64702102

Description:	MOUSE ANTI HUMAN CD80
Specificity:	CD80
Other names:	B7-1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	rA11B11-01/2A7
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant IgG containing a SpyTag at the C-terminus of the antibody heavy chain – liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen Recombinant CD80 (Val35-Asn242) expressed by mammalian cells

External Database

Links

UniProt:

[P33681](#) [Related reagents](#)

Entrez Gene:

[941](#) CD80 [Related reagents](#)

Synonyms

CD28LG, CD28LG1, LAB7

Specificity

Mouse anti Human CD80 antibody, clone rA11B11-01/2A7, is a recombinant monoclonal antibody that recognizes CD80, also known as B-lymphocyte activation antigen B7 or B7-1.

CD80 is a member of the immunoglobulin superfamily, consisting of an extracellular domain bearing two Ig-like domains, a transmembrane sequence, and a short cytoplasmic domain. CD80 is expressed by macrophages, dendritic cells, T cells and activated B-cells. CD80, along with CD86, act as co-stimulatory molecules and are both ligands for CD28 and CTLA-4. Binding to CD28 induces T-cell proliferation, whereas binding to CTLA-4 inhibits T-cell activation. Although CD80 binds to the same receptors as CD86, it displays different characteristics in its avidity and binding kinetics.

CD80 is also expressed on the surface of tumor cells. It binds PD-L1, promoting T cell responses by interfering in the interaction between tumor cells and T cells. CD80 has been proposed as an immunotherapy target.

Flow Cytometry

Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/TZA077GA>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight@488 , DyLight@550 , DyLight@650 , DyLight@680 , DyLight@800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M446225:251027'

Printed on 09 Jan 2026