

## Datasheet: TZC002

**BATCH NUMBER 156817**

<b>Description:</b>	BiSpyCatcher2
<b>Name:</b>	BISPYCATCHER2
<b>Format:</b>	Catcher
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Protein/peptide tag
<b>Product Form</b>	BiSpyCatcher2 purified recombinant protein - liquid
<b>Preparation</b>	Recombinant protein expressed in <i>E. coli</i> and purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.1% ProClin 950
<b>Approx. Protein Concentrations</b>	80 µM (2.3 mg/ml)

**Product Information** **BiSpyCatcher2** (H-BiSpyC2) is a 29 kDa homodimer of SpyCatcher2 ([TZC001](#)) bridged with a flexible linker. Each SpyCatcher2 subunit can form a stable covalent isopeptide bond with a second protein that includes a SpyTag. The reaction occurs spontaneously on mixing, is rapid, irreversible, high-yielding and shows good specificity. BiSpyCatcher2 contains a His6-tag at the N-terminus. To avoid deamidation of asparagine, a N105D mutation was introduced ([Hentrich et al. 2020](#)) for each of the two Catchers.

The reaction is robust at pH 5 to 8, at temperatures from +4°C to +37°C, in various buffer conditions (Ca<sup>2+</sup>/Mg<sup>2+</sup> not needed) and in the presence of detergents. The reaction also occurs inside cells (*in vivo*). Since the speed of the coupling reaction is concentration dependent, it is recommended to use BiSpyCatcher2 undiluted.

All the different formats of BiSpyCatcher2 are compatible with Bio-Rad's recombinant HuCAL® antibodies with a SpyTag2 at the C-terminus of the heavy chain.

[Download BiSpyCatcher coupling protocol](#)

BiSpyCatcher2 (TZC002) is available in two modified versions, BiSpyCatcher2-CYS (TZC002CYS) and BiSpyCatcher-CYS3 (TZC002CYS3), which contain one or three cysteine residues, respectively, for site-specific conjugation by the user.

[View all available Catchers](#)

---

<b>Protein Molecular Weight</b>	Predicted 29018 Da
---------------------------------	--------------------

---

<b>Purity</b>	≥90% determined by SDS-PAGE under reducing conditions and visualized by coomassie blue staining
---------------	---

---

<b>Instructions For Use</b>	<a href="#">View BiSpyCatcher Coupling Protocol</a>
-----------------------------	---

---

<b>Storage</b>	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Acknowledgements</b>	ProClin is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow Sold under license of patent and patent applications, including but not limited to the international PCT application PCT/GB2018/051065, the European patent EP2534484B1, UK patent GB1706430.4, the US patents US9,547,003, US10,247,727, US20200299358, US patent application US16/274908 and corresponding patents HuCAL® is a registered trademark of MorphoSys AG
-------------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20480 available at: <a href="https://www.bio-rad-antibodies.com/SDS/TZC002">https://www.bio-rad-antibodies.com/SDS/TZC002</a> 20480
--------------------------------------	--

---

<b>Licensed Use</b>	For <i>in vitro</i> research purposes. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
---------------------	--

---

## Related Products

### Recommended Useful Reagents

[BiSpyCatcher2:Biotin \(TZC002B\)](#)

[BiSpyCatcher2-CYS \(TZC002CYS\)](#)

[BiSpyCatcher2-CYS3 \(TZC002CYS3\)](#)

[BiSpyCatcher2:HRP \(TZC002P\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M378921:210312'

Printed on 21 Sep 2023

---