

Datasheet: TZC002

BATCH NUMBER 160509

| | |
|----------------------|---------------------|
| Description: | BiSpyCatcher2 |
| Name: | BISPYCATCHER2 |
| Format: | Catcher |
| Product Type: | Recombinant Protein |
| Quantity: | 1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------------|-----|----|----------------|--------------------|
| Immunoassay | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|---|
| Target Species | Protein/peptide tag |
| Product Form | BiSpyCatcher2 purified recombinant protein - liquid |
| Preparation | Recombinant protein expressed in <i>E. coli</i> and purified by affinity chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.1% ProClin 950 |
| Approx. Protein Concentrations | 80 µM (2.3 mg/ml) |

Product Information **BiSpyCatcher2** (H-BiSpyC2) is a 29 kDa homodimer of SpyCatcher2 ([TZC001](#)) bridged with a flexible linker. Each SpyCatcher2 subunit can form a stable covalent isopeptide bond with a second protein that includes a SpyTag. The reaction occurs spontaneously on mixing, is rapid, irreversible, high-yielding and shows good specificity. BiSpyCatcher2 contains a His6-tag at the N-terminus. To avoid deamidation of asparagine, a N105D mutation was introduced ([Hentrich et al. 2020](#)) for each of the two Catchers.

The reaction is robust at pH 5 to 8, at temperatures from +4°C to +37°C, in various buffer conditions (Ca²⁺/Mg²⁺ not needed) and in the presence of detergents. The reaction also occurs inside cells (*in vivo*). Since the speed of the coupling reaction is concentration dependent, it is recommended to use BiSpyCatcher2 undiluted.

All the different formats of BiSpyCatcher2 are compatible with Bio-Rad's recombinant HuCAL® antibodies with a SpyTag2 at the C-terminus of the heavy chain.

[Download BiSpyCatcher coupling protocol](#)

BiSpyCatcher2 (TZC002) is available in two modified versions, BiSpyCatcher2-CYS (TZC002CYS) and BiSpyCatcher-CYS3 (TZC002CYS3), which contain one or three cysteine residues, respectively, for site-specific conjugation by the user.

BiSpyCatcher2-CYS can dimerize by formation of a disulfide bond via the free cysteines. Before conjugation to these cysteines, the BiSpyCatcher2-CYS must be reduced e.g., by addition of 5mM DTT and incubation for 1 hr at room temperature, followed by a fast DTT removal step, e.g., by size exclusion chromatography, as DTT can interfere with the conjugation chemistry.

BiSpyCatcher2-CYS3 contains DTT to avoid oligomerization through disulfide bond formation. Immediately before conjugation, to avoid oxidation and disulfide bond formation, carry out a fast DTT removal step, e.g., by size exclusion chromatography.

[View all available Catchers](#)

| | |
|---------------------------------|--------------------|
| Protein Molecular Weight | Predicted 29018 Da |
|---------------------------------|--------------------|

| | |
|---------------|---|
| Purity | ≥90% determined by SDS-PAGE under reducing conditions and visualized by coomassie blue staining |
|---------------|---|

| | |
|-----------------------------|---|
| Instructions For Use | View BiSpyCatcher Coupling Protocol |
|-----------------------------|---|

| | |
|----------------|---|
| Storage | Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use. |
|----------------|---|

| | |
|------------------|---------------------------------|
| Guarantee | 12 months from date of despatch |
|------------------|---------------------------------|

| | |
|-------------------------|--|
| Acknowledgements | ProClin is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. |
|-------------------------|--|

Health And Safety Information Material Safety Datasheet documentation #20480 available at:
<https://www.bio-rad-antibodies.com/SDS/TZC00220480>

Licensed Use For *in vitro* research purposes. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[BiSpyCatcher2:Biotin \(TZC002B\)](#)

[BiSpyCatcher2-CYS \(TZC002CYS\)](#)

[BiSpyCatcher2-CYS3 \(TZC002CYS3\)](#)

[BiSpyCatcher2:HRP \(TZC002P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391880:211015'

Printed on 21 Sep 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)