

Datasheet: TZC020

BATCH NUMBER 165886

Description:	HUMAN IgA-FcSpyCatcher3
Name:	HUMAN IgA-FcSPYCATCHER3
Format:	Catcher
Product Type:	Recombinant Protein
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Protein/peptide tag
Product Form	FcSpyCatcher3 purified recombinant protein, fused to the hinge region, CH2, and CH3 of human IgA - liquid
Preparation	Recombinant protein expressed in mammalian cell line and purified by affinity chromatography
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.0095% MIT
Approx. Protein Concentrations	80 µM (6.4 mg/ml)
Product Information	Human IgA-FcSpyCatcher3 (hIgA-FcSpyC3) contains a SpyCatcher3 (Keeble et al. 2019) fused to the hinge region, CH2 and CH3 of human IgA (Hentrich et al. 2020),

resulting in a 79.7 kDa protein for the disulfide-bridged dimer. Each SpyCatcher3 subunit can form a stable covalent isopeptide bond with a second protein that includes a SpyTag. The reaction occurs spontaneously on mixing, is rapid, irreversible, high-yielding and shows good specificity. The reaction is robust at pH 5 to 8, at temperatures from +4°C to +37°C, in various buffer conditions (Ca²⁺/Mg²⁺ not needed) and in the presence of detergents. The reaction also occurs inside cells (*in vivo*). Since the speed of the coupling reaction is concentration dependent, it is recommended to use hIgA-FcCatcher3 undiluted.

Human IgA-FcSpyCatcher3 is compatible with Bio-Rad's recombinant HuCAL antibodies with a SpyTag2 at the C-terminus of the heavy chain.

[Download FcSpyCatcher coupling protocol](#)

[View all available Catchers](#)

Protein Molecular Weight	Predicted 79656 Da
---------------------------------	--------------------

Purity	≥90% determined by SDS-PAGE under reducing conditions and visualized by coomassie blue staining
---------------	---

Instructions For Use	View FcSpyCatcher Coupling Protocol
-----------------------------	---

Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
-------------------------	--

Health And Safety Information	Material Safety Datasheet documentation #20480 available at: https://www.bio-rad-antibodies.com/SDS/TZC020 20480
--------------------------------------	---

Licensed Use	For <i>in vitro</i> research purposes. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
---------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[SpyTag3 PEPTIDE \(BLP086\)](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M419886:230622'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)