

Datasheet: TZC025CYS

BATCH NUMBER 156816

Description:	SpyCatcher3-CYS
Name:	SPYCATCHER3
Format:	Catcher-CYS
Product Type:	Recombinant Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Protein/peptide tag
Product Form	SpyCatcher3 purified recombinant protein, with an engineered cysteine residue; used for site-specific chemical conjugation to a label of choice - liquid
Preparation	Recombinant protein expressed in <i>E. coli</i> and purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative	0.1% ProClin 950
Stabilisers	5 mM EDTA
Approx. Protein Concentrations	455 µM (7 mg/ml)

Product Information **SpyCatcher3** (H-SpyC3) is a 15.2 kDa protein that forms a stable covalent isopeptide bond with a second protein that includes a SpyTag ([Keeble et al. 2019](#)). The reaction occurs spontaneously on mixing, is rapid, irreversible, high-yielding and shows good specificity. SpyCatcher3 contains a His6-tag at the N-terminus. The reaction is robust at pH 5 to 8, at temperatures from +4°C to +37°C, in various buffer conditions (Ca²⁺/Mg²⁺

not needed) and in the presence of detergents. The reaction also occurs inside cells (*in vivo*). Since the speed of the coupling reaction is concentration dependent, it is recommended to use SpyCatcher3 undiluted.

All the different formats of SpyCatcher3 are compatible with Bio-Rad's recombinant HuCAL® antibodies with a SpyTag2 at the C-terminus of the heavy chain.

SpyCatcher3 couples faster to SpyTag than the earlier SpyCatcher versions. At neat concentrations the difference in speed between SpyCatcher3 and SpyCatcher2 is negligible but with diluted SpyTag and/or SpyCatcher using SpyCatcher3 shows some benefits. Therefore, Spycatcher3 is recommended for immobilization on plates or beads to bind a capture antibody or other proteins with a SpyTag.

[Download SpyCatcher coupling protocol](#)

SpyCatcher3 (TZC025) is available as a modified version, SpyCatcher3-CYS (TZC025CYS), which contains a cysteine residue for site-specific conjugation by the user

[View all available Catchers](#)

Protein Molecular Weight	Predicted 15368 Da
---------------------------------	--------------------

Purity	≥90% determined by SDS-PAGE under reducing conditions and visualized by coomassie blue staining
---------------	---

Instructions For Use	View SpyCatcher Coupling Protocol
-----------------------------	---

Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	ProClin is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow Sold under license of patent and patent applications, including but not limited to the European patent EP2534484B1, UK patent GB1903479.3, the US patents US9,547,003, US10,247,727, US patent application US16/274908 and corresponding patents HuCAL® is a registered trademark of MorphoSys AG
-------------------------	--

Health And Safety Information	Material Safety Datasheet documentation #20484 available at: https://www.bio-rad-antibodies.com/SDS/TZC025CYS 20484
--------------------------------------	--

Licensed Use	For <i>in vitro</i> research purposes. Any re-sale in any form or any other commercial
---------------------	--

application needs a written agreement with Bio-Rad.

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[SpyCatcher3 \(TZC025\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M378926:210312'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)