

Datasheet: VMA00011

Description:	HUMAN ANTI CYTOKERATIN 18/19
Specificity:	CYTOKERATIN 18/19
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	AbD03748
Isotype:	HuCAL Fab bivalent
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	A bivalent human recombinant Fab selected from the HuCAL GOLD phage display library, expressed in <i>E. coli</i> . This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This product is supplied as a liquid.
Preparation	Recombinant human monoclonal antibody purified by metal chelate affinity chromatography

Source	E.coli
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant human cytokeratin 18
External Database Links	<p>UniProt:</p> <p>P05783 Related reagents</p> <p>P08727 Related reagents</p> <p>Entrez Gene:</p> <p>3875 KRT18 Related reagents</p> <p>3880 KRT19 Related reagents</p>
Synonyms	CYK18
Specificity	<p>Human anti Human cytokeratin 18/19 antibody recognizes human cytokeratin 18, a amino acid type I intermediate filament protein expressed by epithelial cells. The intermediate filaments are expressed as heterodimers in epithelial tissues. Cytokeratin 18 (CK18) is often co-expressed with cytokeratin 8 (CK8). The CK18/CK8 heterodimer is the major cytokeratin pair found in simple epithelia such as that found in the liver, intestine and breast ducts. CK18 plays an important role in maintaining epithelial cell integrity.</p> <p>This clone has also been demonstrated to cross-react with cytokeratin 19. Human anti Human cytokeratin 18/19 antibody recognizes cytokeratin 18 as a band of ~48 kDa and cytokeratin 19 as a band of ~44 kDa in western blotting under reducing conditions.</p>
Western Blotting	Anti cytokeratin 18/19 antibody detects bands of approximately 44 kDa in HeLa cell lysates.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the protein. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>PrecisionAb is a trademark of Bio-Rad Laboratories.</p> <p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See</p>

bio-rad.com/en-us/trademarks for details.
His-tag is a registered trademark of EMD Biosciences.

Health And Safety Information Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/VMA00011>
Antibody (10040)

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Goat Anti Human IgG F(ab')₂ (STAR126...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M429303:240410'

Printed on 10 Apr 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)